



thermo**scientific**

Thermo

ProSightPD

User Guide

Version 4.1 for Proteome Discoverer 2.5

XCALI-98449 Revision A • August 2021

ThermoFisher
SCIENTIFIC

© 2021 Thermo Fisher Scientific Inc. All rights reserved.

Proteome Discoverer is a trademark and Thermo Scientific is a registered trademark of Thermo Fisher Scientific Inc. in the United States. ProSightPC and ProSightPD are registered trademarks of Proteinaceous, Inc. in the United States.

All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

The following are registered trademarks in the United States and other countries:

Windows and Microsoft are registered trademarks of Microsoft Corporation. UniProt is a trademark of European Molecular Biology Laboratory Incorporated Association. TDViewer is a trademark of Tianjin Tandy Digital Technology Co. Ltd. RESID is a registered service mark of John S. Garavelli (individual) in the United States.

Thermo Fisher Scientific Inc. provides this document to its customers with a product purchase to use in the product operation. This document is copyright protected and any reproduction of the whole or any part of this document is strictly prohibited, except with the written authorization of Thermo Fisher Scientific Inc.

The contents of this document are subject to change without notice. All technical information in this document is for reference purposes only. System configurations and specifications in this document supersede all previous information received by the purchaser.

This document is not part of any sales contract between Thermo Fisher Scientific Inc. and a purchaser. This document shall in no way govern or modify any Terms and Conditions of Sale, which Terms and Conditions of Sale shall govern all conflicting information between the two documents.

Revision A, August 2021

Software version: ProSightPD 4.1 for Proteome Discoverer 2.5

For Research Use Only. Not for use in diagnostic procedures.

Contents

	Preface	6
	Access documentation	6
	About this guide	7
	Software installation and activation	7
	Download the software	7
	Activate the software license	8
	Enter the product ID and the activation code	8
	Complete license activation on a computer that is connected to the internet ..	10
	Special notices	11
	Contact us	11
Chapter 1	Overview	12
	Introduction to ProSightPD	12
	Introduction to top-down proteomics	13
	Working with example data	14
	ProSightPD database search algorithms	15
	Annotated proteoform search	16
	Subsequence search	17
	Delta M Mode searches	19
	Label-free quantitation and sliding window deconvolution	20
	ProSightPD task overview	22
	New features in this version	23
Chapter 2	Creating and maintaining databases	25
	Top-down databases in ProSightPD	25
	Create a database in ProSightPD	25
	Download a ready-to-use database	26
	Create a database from a FASTA format file	31
	Change the isoform mass threshold	31
Chapter 3	Working with studies	33
	Studies	33
	Label-free quantitation	33
	Create a study	34
	Process workflows, consensus workflows, and analyses	35

Configure a workflow	36
Add or change an analysis template	36
Add or change a workflow after creating the study.	36
Replace a node in the workflow.	37
Perform Annotated Proteoform searches in the processing workflow.	38
Advanced parameters of the Annotated Proteoform search	40
Perform a Subsequence search in the processing workflow.	41
Advanced parameters of the Subsequence search	43
Perform Single Proteoform searches in the processing workflow	44
Chapter 4 Working with results	46
The results report	46
Proteins page	47
Isoforms page	49
Proteoforms page.	50
PrSMs page	51
MS/MS Spectrum Info page	53
Input Files page	54
Specialized Traces page	54
Quantitative results in ProSightPD	55
Feature Group page.	55
Consensus Feature Group page	56
Quantitative results columns in the Proteoforms page.	57
Visualizing top-down results	59
Use ProSight Lite	59
Use TDViewer	60
Use TDValidator Lite	60
Filter results.	61
Visualize results with charts and graphs	61
Histograms	62
Abundances.	62
Volcano Plot	63
Scatter Plot	64
Chromatogram Traces view.	64
Chapter 5 Workflow templates and analysis templates	66
Chapter 6 ProSightPD nodes	69
Feature Detector nodes	69
Med Res. Feature Detector	69
Hi Res. Feature Detector.	72

Search nodes	79
Annotated Proteoform Search node	79
Single Proteoform Search node	81
Subsequence Search node	83
Tag to Annotated Proteoform Search node	85
Tag to Subsequence Search	88
ProSightPC PUF Writer node	90
cRAWler nodes	90
APD cRAWler node	90
High/High cRAWler node	92
Med/High cRAWler node	95
ProSightPD Consensus nodes	97
PrSM Grouper node	98
Protein Grouper node	98
FDR node	99
Proteoform Validator node	99
Proteoform and Protein Filter node	99
Feature Mapper node	100
Quantifier node	101
PFR Annotator node	101
Chapter 7 Interpreting results	102
P-Score	102
C-Score	103
Expectation Value (E-Value)	103

Preface

This guide describes how to use the Thermo ProSightPD™ 4.1 nodes in the Proteome Discoverer™ 2.5 application to perform top-down proteomics data analysis.

Contents

- [Access documentation](#)
- [Software installation and activation](#)
- [Download the software](#)
- [Activate the software license](#)
- [Special notices](#)
- [Contact us](#)

Access documentation

The ProSightPD application includes complete documentation.

❖ To view user documentation from the Thermo Fisher website

1. Go to Life Sciences Mass Spectrometry Software Download and Licensing Portal <https://thermo.flexnetoperations.com>.
2. Log in to your account.
3. In the navigation pane, click the **Product Search** link.
4. Type **ProSightPD 4.1**.
5. In the Product Search Results, select the entry for *ProSightPD 4.1 User Guide* in the description.

The related software and documentation appear.

6. In the Results column, select the *ProSightPD 4.1 User Guide*.
7. On the Product Download page, select the link under File Name.

8. Click **Download Selected Files**.

About this guide

The *ProSightPD User Guide* describes features specific to top-down analysis. For information about the Proteome Discoverer application, on which the ProSightPD application relies, refer to the *Proteome Discoverer User Guide*.

Software installation and activation

ProSightPD 4.1 is a standalone installer that includes the Proteome Discoverer application and ProSightPD 4.1 nodes. After installing ProSightPD, activate your ProSightPD license in the License Manager tool in the Proteome Discoverer application.

If you have not received your license, contact ThermoMSLicensing@thermofisher.com. Include your sales order number in the email.

Alternatively, you can check your available licenses (and view their activation codes) on the Life Sciences Mass Spectrometry Software Download and Licensing Portal at <https://thermo.flexnetoperations.com>.

Download the software

If you have a license, but do not have your nodes installed, complete the following:

❖ To download the ProSightPD software

1. Go to Life Sciences Mass Spectrometry Software Download and Licensing Portal: <https://thermo.flexnetoperations.com>.
2. Log in to your account.
3. In the navigation pane, click the **Product Search** link.
4. Type **ProSightPD 4.1**.
5. In the Product Search Results, select the entry for *ProSightPD 4.1*.
6. Under File Name (rightside column), click the file name link.

After you download and install the software, activate the license with the License Management interface.

If you have issues with licensing, send an email to ThermoMSLicensing@thermofisher.com. Include a short description of the problem and the ProSightPD license key.

To activate the license, you must have an activation code and product ID number from Thermo Fisher Scientific. Before you transfer a license to another computer, deactivate the license. Refer to the *Proteome Discoverer User Guide* for more information.

Activate the software license

To activate your license, follow these topics as necessary:

1. [Enter the product ID and the activation code](#)
2. [Complete license activation on a computer that is connected to the internet](#)

Enter the product ID and the activation code

To activate your license, you need to know the product ID (XCALI-XXXXX) and the activation code. You can obtain these from the email that Thermo MS Licensing sends you within one week of ordering the software. The email has the subject line “Your Order Is Ready.”

❖ To enter the licensing information

1. Open the Proteome Discoverer application.
2. Choose **Administration > Manage Licenses**.

The Administration page opens.

3. Click **Activate**.

The License Activation dialog box opens to the Activation Code view.

Figure 1. License Activation dialog box opened to the Activation Code view

The screenshot shows a dialog box titled "License Activation" with a close button (X) in the top right corner. The main heading is "Activation Code". Below the heading, there is a paragraph of instructions: "Locate the email from Thermo MS Licensing with 'Your Order is Ready' in the subject line. This email contains the product ID and activation code for your software. If you cannot locate this email, click **Help** for more information." Below this is another paragraph: "Enter the company name, your full name and email address, and the product ID and activation code. Then, do one of the following:" followed by two bullet points: "• If this computer connects to the Internet, click **Online Activation**." and "• If this computer does not connect to the Internet, click **Offline Activation** and follow the instructions. After you retrieve the response file (activation.xml), click **Process Response File** on this computer to complete the activation process." Below the instructions are five input fields: "Company:", "Full Name:", "User Email:", "Product ID:" (with "XCALI-____" pre-filled), and "Activation Code:" (with "____-____-____-____" pre-filled). At the bottom of the dialog box are five buttons: "Help", "Online Activation", "Offline Activation", "Process Response File", and "Cancel".

4. If you have not received your activation code, do the following:
 - a. Check your Junk Email folder.

If the email is not in your Junk Email folder, log in to your account at the following URL: <https://thermo.flexnetoperations.com>.

In the left navigation pane, under Software & Services, click **Order History**. Then, in the list of ordered products, click the order number.

- b. If you cannot find your account, send an email message to Licensing at ThermoMSLicensing@thermofisher.com.

Provide the following information in the body of the message:

- Software application: ProSightPD
- Sales order number or purchase order number
- End user name
- End user email

5. In the License Activation dialog box (Figure 1), enter the following:

- Your company name.
- Your full name.
- Your contact email address.
- The product ID for the ProSightPD application.
- The activation code. You can type or paste the activation code.

Table 1. Product IDs for the ProSightPD application

Material Order No.	Product ID	Description
OPTON-31081	XCALI-65179	SW, PROSIGHTPD 4.1
OPTON-31082	XCALI-65180	SW, PROSIGHTPD 3.0 TO PROSIGHTPD 4.1
OPTON-31083	XCALI-65181	SW, PROSIGHTPD BEFORE 3.0 TO PROSIGHTPD 4.1
OPTON-31084	XCALI-65182	SW, PROSIGHTPD 4.1 UPG FRM PROSIGHTPC

6. Continue to the following topic:

- [Complete license activation on a computer that is connected to the internet](#)

Complete license activation on a computer that is connected to the internet

Follow these instructions if your processing computer has an Internet connection.

❖ To complete license activation on an online computer

1. If you have not already entered the licensing information, enter it in the Activation Code view of the License Activation dialog box.
2. Click **Online Activation** to process the activation code.
3. To close the dialog box, click **OK**.

This completes the online license activation process.

For additional information about deactivating a license or using an offline computer, refer to the *Proteome Discoverer User Guide*.

Special notices


Make sure you follow the special notices presented in this guide. Special notices appear in boxes; those concerning safety or possible system damage also have corresponding caution symbols.

IMPORTANT Highlights information necessary to prevent damage to software, loss of data, or invalid test results; or might contain information that is critical for optimal performance of the system.

Note Highlights information of general interest.

Tip Highlights helpful information that can make a task easier.

Contact us

Contact	Email	Telephone	QR Code
ProSightPD Support	prosightpdsupport@thermofisher.com		
Proteinaceous	info@proteinaceous.net		
U.S. Customer Service and Sales	us.customer-support.analyze@thermofisher.com	(U.S.) 1 (800) 532-4752	

Overview

This chapter introduces the ProSightPD application and provides a brief discussion of top-down proteomics.

Contents

- [Introduction to ProSightPD](#)
- [Working with example data](#)
- [ProSightPD database search algorithms](#)
- [Label-free quantitation and sliding window deconvolution](#)
- [ProSightPD task overview](#)
- [New features in this version](#)

Introduction to ProSightPD

The ProSightPD application is a set of nodes in the Proteome Discover platform. These nodes are designed to analyze top-down mass spectrometry data.

The ProSightPD nodes identify, quantify, and characterize proteoforms. The ProSightPD application is fully integrated into the Proteome Discoverer framework and utilizes a common user interface and results architecture. The primary difference between the ProSightPD nodes and the Proteome Discoverer nodes is the focus on proteoform identification and not peptide identification. Many of the notable differences such as database structure, tables in the results, and graphical result views stem from this fundamental practice.

The ProSightPD application supports several top-down workflows including:

- High-throughput LC/MS proteoform discovery
- Large proteoform discovery using proton transfer charge reduction (PTCR)
- Label-free quantitation (LFQ) of proteoforms
- Targeted proteoform searches


- Analysis of infused samples and manually collected data (not using a predefined MS method)
- Analysis of isotopically and non-isotopically resolved precursors

Introduction to top-down proteomics

The application works with the mass values inferred from mass spectral data from middle-down/bottom-up and top-down proteomics MS/MS experiments.

Figure 2 compares bottom-up and top-down methods.

Figure 2. Comparing top-down and bottom-up proteomics



In top-down proteomics, proteins are introduced into a mass analyzer, where they are subjected to gas-phase fragmentation. The right side of Figure 2 illustrates this process.

Top-down proteomics is a technique for protein identification and characterization. Combining top-down proteomics searches with the shotgun annotation process provides a method for rapid and accurate protein definition.

Top-down proteomics can provide thorough sequence coverage of proteoforms, supporting characterization of complex combinations of PTMs.

Separation and ionization of intact proteins present many challenges beyond the scope of this guide. The bioinformatics of top-down proteomics, though, has some challenges that the application addresses.

- Because the precursor ions are large, they are almost always multiply charged. This complicates spectral comparison techniques used in certain other proteomic strategies. Fortunately, this issue is avoided when converting to neutral masses. After, you collect MS and MS/MS spectra, you sum the relevant scans and then run through an automated analysis to infer mass, using the resulting mass values for protein identification and characterization.
- Because the precursor mass can represent either a highly modified protein or an internal fragment of the intact protein, no single strategy of comparing the observed mass values to a proteome database is guaranteed to identify the protein. For this reason, the application provides the search modes described in [ProSightPD database search algorithms](#).

Working with example data

The ProSightPD application installs with a database and analysis workflow template for the Thermo Scientific™ Pierce™ Intact Protein Standard Mix. You can download an example of the Pierce Intact Protein Standard Mix data file at the Life Sciences Mass Spectrometry Software Download and Licensing Portal (<https://thermo.flexnetoperations.com>).

The Pierce Intact Protein Standard Mix contains six recombinant proteins. For additional information, refer to

<https://www.thermofisher.com/order/catalog/product/A33527#/A33527>.

This data file was collected on a Thermo Fisher Scientific Fusion Lumos Tribrid Mass Spectrometer with high resolution precursor and HCD fragment ions scans.

To become familiar with the ProSightPD data analysis process using example data, complete the following:

1. Download the `Pierce_Intact_Protein_Standard_Top_Down_Example.raw` data file at the Life Sciences Mass Spectrometry Software Download and Licensing Portal (<https://thermo.flexnetoperations.com>).
2. Open the Proteome Discoverer application.
3. On the Start page, click **New Study/Analysis**.

The New Study and Analysis dialog box opens.

4. Type the Study Name.
5. For the Study Root Directory, click the browse icon and locate the path.
6. Click **OK**.

The study opens.

7. Click **Open Analysis Template** and locate the PSPD HiHi TopDownStandard Analysis Template (C:\Users\Public\Documents\Thermo\Proteome Discoverer 2.5\Common Templates\ProSightPD 4.1\Analysis Templates).
8. Click **Open**.
9. Click the **Workflows** tab on the Study page.
10. In the Processing Step area of the Analysis pane, click **Edit**.
11. In the Workflow Tree, select the Annotated Proteoform Search node.
12. In the Parameters of the 'ProSightPD 4.1 Annotated Proteoform Search' area on the left side of the Workflow Tree, click the Input Database dropdown arrow and select the **PierceIntactProteinStandardMix** database file (it is the only database available for newly installed software).
13. Click the **Input Files** tab on the Study page.
14. Click **Add Files**.

The Add Files dialog box opens.
15. Locate the `Pierce_Intact_Protein_Standard_Top_Down_Example.raw` data file that you downloaded at the beginning of this process, then click **Open**.

The file is now listed on the Input Files page.
16. Drag and drop the file into the Processing Step area of the Analysis pane on the right side.

The Run button now activates and turns green.
17. Click **Run**.

Once you click Run, the Administration page automatically opens where you can monitor the progress of your analysis. For additional information about viewing your results, see [Chapter 4, "Working with results."](#)

ProSightPD database search algorithms

The following topics describe the search algorithms specific to the ProSightPD application. Each search mode overcomes different issues of protein identification and characterization:

- [Annotated proteoform search](#)
- [Subsequence search](#)
- [Delta M Mode searches](#)


Annotated proteoform search

The annotated proteoforms search matches MS/MS data against all intact proteoforms of proteins in a database. It is the defining search mode for top-down proteomics. Annotated proteoform searches use the precursor mass to generate a subset of the proteome database to query.


For each proteoform with a theoretical precursor mass within the window of the observed precursor ion mass, plus or minus the defined search tolerance, the annotated proteoform search compares all theoretical fragments and masses to observed fragment ion masses. The application determines the number of observed fragment ions matching the fragment tolerance and uses this value to score the identification (Figure 3).

The application calculates scores for each probable MS/MS to candidate proteoform sequence match (PrSM). An observed ion matches a theoretical ion when the two masses are within a user-defined tolerance. The probability of the observed number of fragment ions matching by chance is then determined and reported as a P-Score.

Figure 3. Annotated proteoforms searches



Although the application queries each proteoform with a theoretical precursor mass in the window, it displays only those proteoforms that meet the user-defined filtering of search results. These filters include the minimum number of matching fragments and false discovery rate (FDR) (Figure 4).

Figure 4. Annotated proteoforms search process

Use the following strategies when running annotated proteoform searches:

- Use a 1000-Da precursor search window search as the first search for an unknown protein. If there are few modifications on the unknown protein not in the Input Database, a 1000-Da intact search frequently identifies, but does not characterize, the protein.

A large number of ions matching one terminal in a protein is evidence of a protein's identity. You can use the ProSight Lite and the TDValidator Lite applications for further characterization. For more information on these tools, see [Use ProSight Lite](#) and [Use TDValidator Lite](#).

- Use Delta M (Δm) Mode to locate unknown modifications near either terminus. If the 1000-Da annotated proteoform search fails to identify a protein, consider running another annotated proteoform search with a 1000-Da precursor search window with Delta M Mode enabled.

Activating Delta M Mode increases the likelihood that the search will identify proteins with unknown modifications. However, this mode takes approximately two times longer than the corresponding annotated proteoform search. For more information, see [Delta M Mode searches](#).

Subsequence search

A ProSightPD subsequence search matches MS/MS data against all subsequences of all proteoforms of proteins in a database. The subsequence search is similar to a bottom-up no enzyme search. A subsequence search is a “brute force” search of an entire database and can take a long time. It looks at every possible subsequence of every base proteoform (unless

mentioned otherwise) in the database and attempts to identify any subsequence that matches the observed intact ion mass within a tolerance. For each subsequence matching the intact ion mass, the subsequence search performs an absolute mass search and reports any subsequence that matches the observed intact ion mass and is able to generate the observed fragment ion pattern.


In a typical top-down experiment, not all of the polypeptides identified are intact proteins. A subsequence search identifies those proteins that are a product of biological degradation and cannot be logically predicted. This search compares the observed precursor mass to all possible entries of a particular database within a defined tolerance, for example, less than 10 ppm. The application theoretically fragments those entries that fall within the defined tolerance and compares the observed fragment ions.

- Use a subsequence search if an absolute mass analysis fails to identify a protein.
- The default subsequence search searches only for the basic proteoforms with no known modifications. However, searching for modified proteoforms increases the search run time.

A subsequence search is a two-step process that is repeated for each base protein sequence in the proteome database:

1. Identify a candidate entry matching an observed precursor mass.
2. Calculate all possible theoretical fragment ions for the candidate entry, and then compare the theoretical fragment ion masses to the observed fragment ion masses.

Figure 5 shows the process involved in a subsequence search.

Figure 5. Subsequence searches

In a subsequence search, the precursor search tolerance is an estimate of measurement error on the observed precursor mass. The value is usually small compared to an absolute mass precursor search window.

Delta M Mode searches

Delta M (Δm) Mode is a feature for identifying proteoforms containing unknown mass shifts, for example, PTMs. The delta is the difference between the observed precursor mass and the theoretical precursor mass. When you perform a search in Delta M Mode, the ProSightPD software concurrently performs three queries per sequence to compare the following:

- The theoretical fragment ion masses of the protein sequence to the observed fragment ion list as usual
- The theoretical fragment ion masses derived from the sequence and the Delta M applied N terminal to the observed fragment ion mass list
- The theoretical fragment ion masses derived from the sequence and the Delta M applied C terminal to the observed fragment mass list


A Delta M Mode search takes approximately two times longer than the same search without Delta M Mode.

By carefully observing the pattern of fragments with and without Delta M (Δm) Mode, you can frequently locate the source of the mass shift.

For example, a particular result returns with the two smallest N-terminal fragments matching without the mass shift. If all other matching N-terminal fragments contain the mass shift, the unknown mass can be localized on an amino acid between the second and third N-terminal fragments (Figure 6).

When Delta M Mode is enabled, fragments found using the Delta M Mode display on the PrSMs page, along with two additional columns: Delta M Residue and Delta M Fragment(s), where the delta mass and the possible residue where the delta mass was localized. Delta masses also display in the ProSightPD fragment map with an orange square around the potential residue where the delta mass was localized.

Figure 6. Schematic of Δm mode



Label-free quantitation and sliding window deconvolution

The ProSightPD application performs label-free quantitation (LFQ) using the sliding window deconvolution algorithm.


In an LFQ experiment, the proteoforms are quantified by elution profile, called Feature Groups. A Feature Group is the decharged and deisotoped signal measured over time in an LC/MS experiment. Feature Groups are measured in the Feature detector nodes using the sliding window algorithm in conjunction with the Xtract (for high resolution, isotopically-resolved data) or kDecon (for low or medium resolution, isotopically unresolved data) algorithms for deconvolution.

The sliding window algorithm averages spectra over a succession of windows in retention time, deconvolves each average spectrum, and then merges similar masses from consecutive deconvolutions to form feature groups.

Figure 7 illustrates the concept of sliding windows.

Note FWHM refers to full-width half-max.

Figure 7. Sliding window overview



The sliding window algorithm benefits top-down LFQ in several ways:

- Reduces the number of false positives due to noisy data
- Improves sensitivity using signal averaging
- Identifies co-eluting species better
- Defines elution profiles for quantitation

The ProSightPD Hi Res. Feature Detector node and the ProSightPD Med. Res. Feature Detector node contain the parameters controlling the sliding window algorithm.


1. After the detector node determines the feature groups for data files, the node maps the feature groups to individual PrSMs and connects feature groups to quantitation traces. The node groups the features that are similar in mass and retention time between files. A group of feature groups across all files is a *consensus feature group*.
2. The Feature Mapper node maps the consensus feature groups to proteoforms based on a mass and retention time threshold.

3. The ProSightPD Quantifier node produces quantification ratios and statistics using the measured abundances and study factors you applied.

ProSightPD task overview

Follow the instructions in [Figure 8](#) to specify preferences and perform the search.

Figure 8. The ProSightPD task flow



New features in this version

The new and enhanced features in this release include:

- Custom Mass Modifications
 - Ability to add a Custom Mass Modification
 - Ability to enter user-defined Monoisotopic mass
- Locking feature for PTMs
 - Greater control over database size
 - Improved search speeds
 - Ensures all search candidates include the locked modifications at search time
- Delta M Mode reporting
 - Reports delta mass
 - Reports delta mass fragments and potential location of delta mass
 - Enhanced delta m mode visualization
- MzIdentML into Database Manager
 - Proteome Discoverer Results in MzIdentML are compatible inputs for Database Manager
 - Use bottom-up experimental results to create TopDown databases
- Heavy Labeled search option
 - Support for C13 and N15 labeled intact proteins
 - Utilizes custom averaging tables
 - Assumes complete labeling
- Database dropdown menu
 - You can now select a database from a dropdown menu in the search node.

Note Only databases in the ProSightPD database folder are available in the dropdown menu. If a new database is created, you must close and reopen the study to access the new database.

- MassShift/Recal in cRAWler
 - Shifts masses by a constant amount in parts per million.
 - The new parameters are in the cRAWler node, and separate mass shifts can be set for the precursor and fragments.
 - This feature accounts for calibration errors, which are constant throughout the data file, and does not re-calibrate the data but rather applies a constant user-defined shift to all masses in all spectra in the data file.

- Labile modification mass shift
 - Lets you input known labile modification masses to be considered during the search.
 - Shifted precursor and fragment masses will be considered during the search.
- MS2/MS3 cRAWler
 - MS3 data can now be analyzed in ProSightPD.
 - To pass MS3 spectra to the cRAWler node for deconvolution and grouping, a new grouper mode has been added: “MS2 & MS3”

Note In the Spectrum Selector node, the Precursor Selection must be set to “Use MS(n-1) Precursor” and the MS Order must include MS2 and MS3.

- PFR accession
 - Each proteoform result receives a Proteoform Reference accession number (PFR). This lets you track proteoforms across experiments.
- Disulfide bond toggle
 - To overcome incorrect disulfide bond annotations.

Creating and maintaining databases

This chapter describes working with top-down databases in the ProSightPD application.

Contents

- [Top-down databases in ProSightPD](#)
- [Create a database in ProSightPD](#)

Top-down databases in ProSightPD

The key input for top-down data analysis is a properly annotated database. Protein databases are available in many formats and from several sources. The ProSightPD application supports only UniProt-formatted extensible markup language (XML) and UniProt-formatted FAST-All (FASTA) file formats. Thermo Fisher Scientific strongly recommends using UniProt-formatted XML files when creating a database as these files include the full complement of known modifications.

You can download protein databases and full proteome databases from:

- UniProt—for a comprehensive source <https://www.uniprot.org>.
Select an XML file.
- Proteinaceous— for selected, common proteome databases in XML and ProSight DataBase (PSDB) formats <https://proteinaceous.squarespace.com/prosightpdfaqs>.

In the ProSightPD software, the final usable database format is PSDB. The Database Manager can convert the XML and FASTA files to PSDB format.


Create a database in ProSightPD

Proteinaceous offers several ready-to-use databases for common organisms, which are in XML format that you can import through the database management tool in the ProSightPD application.

Download a ready-to-use database

Tip If you want only some proteins from a genome or proteins from different genomes, download from UniProt. If you want a complete genome, download from Proteinaceous.

❖ Download and configure an XML format database

- Go to one of the following websites:
 - www.uniprot.org
 - <https://www.proteinaceous.net/>
- Locate and download the XML file.
- Save the file.
- In the Proteome Discoverer application, click the ProSightPD Database Manager icon, , in the toolbar.


Tip You can also access the Database Manager from the Help menu.

The ProSightPD Database Manager dialog box opens.

- In the dialog box, do the following:
 - Click **Create ProSightPD Database**.
 - Click **Browse**, and locate the saved XML file.
 - Select the file.

The database file appears in the dialog box (Figure 9). The application displays the number of proteins found in the imported file.

Figure 9. ProSightPD Database Manager dialog box



If the imported database is correctly formatted, a green check and the number of correctly formatted proteins appear. If the file is incorrectly formatted, a red X appears.

If the number of entries shown is fewer than you expected, the lower number might be due to the Max Isoform Mass setting being too low. For information on changing that setting, see [Change the isoform mass threshold](#).


Other reasons that the number of isoforms might be unexpectedly low:

- A sequence contains an unrecognized amino acid (B,Y,Z, and so forth.) or a wildcard.
- The entry header is incorrectly formatted (that is, not UniProt format).
- Some of the entries were duplicated.

The Database Manager removes incorrectly formatted and duplicate entries.


❖ (Optional) Edit isoforms in Protein Annotator

Figure 10. Protein Annotator view




In the Protein Annotator dialog box, you can do the following:

- Add or remove isoforms from a database.
 - Add or remove modifications from an isoform.
1. To edit or annotate the database, click **Protein Annotator from the ProSightPD Database Manager**.
 2. To modify an isoform:
 - a. Select the isoform.

- b. Select the residue to modify.
- c. Select the desired modification in the Point Features pane. You can lock a feature by selecting the lock icon  next to the modification.


Note After a modification is selected and locked, you cannot add additional modifications to the selected residue. The selected modification will be on all search candidates.

3. To add a disulfide bond or endogenous cleavage:
 - a. Click the plus icon  in the Range Features pane.
The Add Range Features dialog box opens.
 - b. In the Add Range Feature dialog box, select the UniProt Feature Type and type the Modification Start Point and End Point.
 - c. Click the **Save and Exit** icon in the lower right corner.
4. To add a mutation or cSNP (single-nucleotide polymorphism), select the residue to mutate.
 - a. Select the cSNPs under Point Features.
 - b. Select the mutation.
 - c. (Optional) In the protein sequence pane, you can add one or more modifications in the Point Features pane.

Tip To remove a modification, press CTRL and select the modification. To select multiple point features, press the CTRL key.

5. To add Custom Mass Modifications, select the residue to modify.

Figure 11. Custom Mass Modifications







- Choose the Custom Mass Modifications tab in the Point Features pane.
- Type the Custom Mass in Daltons (Da).
- Click the green plus icon to add the modification to the modifications list and apply the modification to the residue. If there is more than one Custom Mass Modification listed, you can select multiple modifications to add to each residue.

Note Unlike named PTMs, the addition of multiple Custom Mass Modifications will result in orange colored squares.

6. Click the **Save and Exit** icon () in the lower right corner.

The downloaded database now appears in the list of searchable databases ([Figure 13](#)).

Figure 12. The ProSightPD Database Manager


Input	Name	Description	# Entries	# Isoforms	Size (MB)	Date
	Yeast.psdB	yeast	6,049	6,077	5.78	06/12/2020
	Xenograft_IP.psdB	IP	990	1,362	1.38	10/12/2020
	UpdatedMouse08182020.psdB	does have doub	13,027	17,544	18.25	08/18/2020
	TopDown Standard.psdB	TopDown Stand	5	5	.04	06/12/2020

❖ **To add an available feature to the included feature set**

1. Click **Modify Available PTMs** .


The Global Feature Set Editor opens (Figure 13).

Figure 13. Global Feature Set Editor



2. If the modification is not available in Included Features, search for it under Available Features.
3. You can apply a modification on a group or an instance level:
 - To make the modification group available, in the Available Features list, click the up arrow for the row to move the complete modification group (that is, all the individual instances) to the Included Features list.
 - To make an instance available, expand the modification group. Select the up arrow on the instance row to move the instance to the Included Features list.

You can add additional instances one at a time.

4. (Optional) To globally apply the modification: in the Included Features list, select the modification group or instance, and click the plus icon ().
5. Click **Save**.

The new modifications are available.

Create a database from a FASTA format file

While XML files are strongly encouraged for database creation, you can also use FASTA files in ProSightPD.

Note Be aware that FASTA files do not contain any PTM information. As such, only N-terminal modifications defined in the search nodes or modifications manually added through Protein Annotator are considered if FASTA files are used.

Make sure the FASTA file format matches the UniProt format, particularly in the format for the file description line. Otherwise, ProSightPD might not be able to process the file.

To enable the software to process the FASTA file, edit the description line to include:

- >sp or >tr
- unique accession for each entry
- pipe (|) before and after the unique accession and no spaces

This is the UniProt format:

```
>sp|P02144|MYG_HUMAN Myogl obi n OS=Homo sapi ens GN=MB PE=1 SV=2
MGLSDGEWQLVLNVWGKVEADI PGHGQEVLI RLFKGGHPETLEKFDKFKHLKSEDEMKASE
DLKKHGATVLTALGGI LKKKGHHEAEI KPLAQSHATKHKI PVKYLEFI SECI I QVLQSKH
PGDFGADAQGAMNKALELFRKDMASNYKELGFQG
```

This is a non-conforming example:

```
>MYG_HUMAN
MGLSDGEWQLVLNVWGKVEADI PGHGQEVLI RLFKGGHPETLEKFDKFKHLKSEDEMKASE
DLKKHGATVLTALGGI LKKKGHHEAEI KPLAQSHATKHKI PVKYLEFI SECI I QVLQSKH
PGDFGADAQGAMNKALELFRKDMASNYKELGFQG
```

This is an example of an acceptable user-generated FASTA file:

```
>sp|P02144|MYG_HUMAN
MGLSDGEWQLVLNVWGKVEADI PGHGQEVLI RLFKGGHPETLEKFDKFKHLKSEDEMKASE
DLKKHGATVLTALGGI LKKKGHHEAEI KPLAQSHATKHKI PVKYLEFI SECI I QVLQSKH
PGDFGADAQGAMNKALELFRKDMASNYKELGFQG
```

Change the isoform mass threshold

This setting limits the number of isoforms included in your database. You can use the isoform mass threshold to improve search speed by limiting the database entries to only isoforms likely to be found in the sample.

For example, a typical top-down LC/MS experiment is unlikely to detect an intact form of Titin (over 1Mda).

Alternatively, the entire yeast proteome includes approximately 6050 entries. However, only 4550 are below the default 70000 threshold.

❖ To change the upper mass limit of isoforms


1. Click the gear icon (⚙️) in the Database Manager.

The Advanced Parameters dialog box opens (Figure 14).

2. For Max Isoform Mass, type a value.

The default is 70 000 Da.

Figure 14. Advanced Parameters dialog box



3. Click the **Save** icon.
4. Continue creating your database. See [Create a database from a FASTA format file](#).

Working with studies

This chapter describes studies and analyses.

Contents

- [Studies](#)
- [Create a study](#)
- [Configure a workflow](#)
- [Advanced parameters of the Annotated Proteoform search](#)
- [Advanced parameters of the Subsequence search](#)

Studies

A study is a container for your data files and the analyses performed with those data files. The study is also where you structure and annotate your data with meta information, such as the biological condition (for example, disease versus control), analysis time points, biological replicate number, and so on.

The ProSightPD application analyzes data by running a processing workflow and a consensus workflow in the study.


Label-free quantitation

The ProSightPD application performs the label-free quantification (LFQ) of proteoforms using these nodes:

- ProSightPD Hi Res Feature Detector or the ProSightPD Med Res Feature Detector node in the processing workflow
- Feature Mapper node and ProSightPD Quantifier node in the consensus workflow

[Figure 15](#) shows the workflows.

Figure 15. Processing workflow and consensus workflow for label-free quantitation



Create a study

The most convenient way to create a new study is to use the New Study and Analysis dialog box. You can add files to the study after creating it, but the import function is available only in the New Study and Analysis dialog box.

❖ To create a study using the New Study and Analysis dialog box

1. On the Start page, click **New Study/Analysis**.

The New Study and Analysis dialog box opens.

Figure 16. New Study and Analysis dialog box

The screenshot shows the 'New Study and Analysis' dialog box with the following fields and values:

- Study Name: My Top Down Study Name
- Study Root Directory: C:\Data
- Import From File: (empty)
- Processing Workflow: ProSightPD 4.0 \ Workflows \ PSPD LFQ for Med HI data.pdProcessingWF
- Consensus Workflow: ProSightPD 4.0 \ Workflows \ PSPD 1 percent FDR Consensus.pdConsens

Buttons for 'OK' and 'Cancel' are visible at the bottom right.

2. Type the Study Name.
3. For the Study Root Directory, click the browse icon and locate the path.

Put the study directory on a local drive that has approximately the same free space as your raw data file size.

4. (Optional) For Import from File, click the browse icon and locate the file that contains the study definitions.
5. For the Processing Workflow, click the browse icon and locate the processing workflow template that matches the purpose of your study.

To use an analysis template, which includes the processing workflow and the consensus workflow, leave the workflow fields empty. The analysis templates are not available until the study is created.

For more information on workflows and templates, see [Process workflows, consensus workflows, and analyses](#).

6. For the Consensus Workflow, click the browse icon and locate the consensus workflow template based on how you want to filter your results.

For more information on workflows and templates, see [Process workflows, consensus workflows, and analyses](#).

Process workflows, consensus workflows, and analyses

A workflow is composed of a series of nodes, and each node performs an individual step of the data analysis. Depending on the workflow's purpose, they usually consist of 5 to 10 nodes. Most nodes have their own parameters to set. The basic functions of the workflows are:

- In the processing workflow, the software searches the database and assigns PrSMs. If the analysis is quantitative, the processing workflow also performs feature detection.
- The consensus workflow filters the preliminary results (for example, by the FDR node) and groups them into proteoforms, isoforms, and proteins. The results display these groups. If there is quantitative analysis in the processing workflow, then the consensus workflow performs the required feature mapping and statistical analysis.

You can run multiple processing workflows separately and then combine them in the consensus workflow for unified reporting.

To simplify creating a workflow, the software includes a set of fit-for-purpose templates. Build your workflow from the provided templates, which have been tested and optimized for confident proteome analysis. Many of the node parameters already have their appropriate settings in the template.

There are processing workflow templates, consensus workflow templates, and analysis templates, which contain both the processing and consensus workflows. The templates are named for their intended use. For example, the Comprehensive Discovery template includes nodes that cover the largest search space.

Configure a workflow

The following sections describe how to add or change workflows and other related tasks:

- [Add or change an analysis template](#)
- [Add or change a workflow after creating the study](#)
- [Replace a node in the workflow](#)
- [Perform Single Proteoform searches in the processing workflow](#)

Add or change an analysis template

❖ To add or change an analysis template

1. Create a study in the dialog box without selecting either a processing or a consensus workflow.

For instructions, see [To create a study using the New Study and Analysis dialog box](#).

2. On the Study Definition page, click **Open Analysis Template** .
3. In the File Explorer, open the analysis template:

C:\Users\Public\Documents\Thermo\Proteome Discoverer *version*\Common Templates\ProSightPD *version*\Analysis Templates

4. Add files, study factors, and a database and adjust the parameters to suit your analysis.

Add or change a workflow after creating the study

Tip You can overwrite the existing workflow with the new workflow, that is, you do not need to delete the existing workflow first.

❖ To add or change a workflow after creating the study

1. Do one of the following:
 - To use a workflow supplied in the application, click **Open Common**.
 - To use a workflow saved on your computer, click **Open**.
2. Locate the workflow, and click **Open**.
3. (Optional) Repeat to add or change another workflow.
4. Add files, study factors, and adjust the node parameters to suit your analysis.

Replace a node in the workflow

You might want to change a workflow node for many reasons, for example, you are getting unexpected results, too many results, or changing from a database to a single proteoform search.

Note The ProSightPD 4.1 application includes the search nodes from the ProSightPD 3.0 application for you to more conveniently reproduce previous results in the current version. These legacy nodes require a database in the .pscw file format generated in ProSightPC. You cannot directly convert .pscw to .psdb. Use the source .xml file, which the .PSCW file was generated from, to generate a .psdb format database. When using these nodes, replace the 4.1 node with the equivalent legacy node and review all parameters to ensure agreement between your ProSightPD 3.0 workflows and your ProSightPD 4.1 workflows.

❖ To replace a node in the workflow

1. Select the node that you want to replace.

This example shows the ProSightPD High/High cRAWler node selected.

The screenshot displays the ProSightPD software interface. On the left, the 'Workflow Nodes' panel lists various search and processing nodes. The 'ProSightPD Search' category is expanded, and the 'ProSightPD High/High cRAWler' node is highlighted with a red dashed box. The main 'Workflow Tree (Reprocess)' pane shows a flowchart with five nodes: 'Spectrum Files' (0), 'Spectrum Selector' (1), 'ProSightPD High/High cRAWler' (2), 'ProSightPD 4.1 Annotated Proteoform' (3), and 'ProSightPD 4.1 Subsequence Search' (5). The 'ProSightPD High/High cRAWler' node is selected, and its connections to the subsequent nodes are visible. The top of the interface shows tabs for 'Study Definition', 'Input Files', 'Samples', 'Analysis Results', 'Workflows', and 'Grouping & Quantification'. The 'Workflows' tab is active, showing the current workflow name and a description field.

2. Press DELETE on the keyboard.
3. Drag the new node into the workflow pane, and connect the nodes.

For example, to analyze top-down data that contains unresolved precursor signals, use the ProSightPD Med Hi cRAWler node.

Perform Annotated Proteoform searches in the processing workflow

The Annotated Proteoform Search node is designed to quickly analyze raw data. This node considers only proteoforms included in your database. This includes any known PTMs, truncations, or amino acid variants that are annotated in the database.

Annotated Proteoform searches can be applied in the following two ways:

- A narrow tolerance search
- A wide tolerance search

A narrow tolerance search is usually conducted with a 2.2 Da precursor mass tolerance to account for the common “off-by-one” data deconvolution artifacts. The narrow tolerance search completes quickly and the results do not require further analysis to improve confidence (a few exceptions might include: loss of hydrogen from S-S bonds and deamidation).

A wide tolerance search typically utilizes a 100-1000 Da precursor mass tolerance. This search returns many proteoforms whose precursor mass varies significantly from the theoretical precursor mass. The goal of this search is to find proteoforms that might have an unknown modification resulting in a larger precursor mass error, which could only be found using a wide tolerance. These results require additional work to determine the modification.

Thermo Fisher Scientific recommends opening individual results in the ProSight Lite application to further explore potential modifications. If the modification is identified, it can be added to the database. Then the search must be re-processed to find all proteoforms containing that modification.

To set advanced parameters for an Annotated Proteoforms search, see [“Advanced parameters of the Annotated Proteoform search.”](#)

❖ To add an Annotated Proteoform search in the processing workflow

1. Open a workflow template that includes the ProSightPD Annotated Proteoform Search node.
 - a. If the workflow does not include the Annotated Proteoform Search node, drag and drop the Annotated Proteoform Search node into the workflow.
 - b. Connect the node by dragging the arrows from the upstream node to the Annotated Proteoform node.
2. Select the Annotated Proteoform node to view the parameters.
3. Set the Input Database by selecting from the dropdown menu.

Note To create a new database launch the Database Manager. To select a newly created database, you must close and reopen any open studies.

4. Set the Precursor Mass Tolerance.

The following are recommended values for narrow or wide tolerance searches:

- Narrow tolerance search is 2.2 Da
- Wide-tolerance search 200 Da

Note If you are think there might be a large modification like a signal peptide, consider setting your tolerance larger than the potential modification.

5. Set the Fragment Mass tolerance.

Thermo Fisher Scientific recommends 10 ppm. Use the narrowest tolerance possible for your data. As the tolerance widens, the search time and false positive rate increases.

6. Set the FDR calculation.

The default is True, which means that the application calculates the FDR based on the threshold set in the FDR node in the consensus workflow. To calculate the FDR, the FDR node in the consensus workflow is required. If your sample is not complex or if you want to see all PrSMs, set this value to False to not calculate an FDR.

Note Results are automatically displayed at the specified FDR cutoff for each result level (PrSM level, Proteoform level, Isoform level, and so on)

7. (Optional) Select a static modification to apply to a specific residue.

Currently, you can only apply one static modification in the search node. To set additional modifications, use the Protein Annotator in the Database Manager.

Advanced parameters of the Annotated Proteoform search

The default parameters of the Annotated Proteoform search are designed to work in most situations, however, most analyses can be improved/customized by optimizing the advanced parameters listed below. For a full list of advanced parameters, see [“Annotated Proteoform Search node.”](#) Table 2 is a list of advanced parameters and a brief description of how they impact the Annotated Proteoform search and results.

Table 2. Advanced parameters for the Annotated Proteoform

Parameter	Description
Maximum PTMs per Isoform	Limits the number of PTMs that are expanded and searched during search time. For example, if an isoform has 10 PTMs in the database and the maximum is set to 4, then only proteoforms with up to 4 PTMs are searched. This improves the search time. However, in cases of highly modified proteoforms, you should increase this value.
Maximum SNPs per Isoform	Limits the number of amino acid mutations of single nucleotide polymorphisms (SNPs) in the same way as the Maximum PTMs per isoform parameter.
Remove Disulfide Bonds	Thermo Fisher Scientific recommends setting this parameter to False. You can set this parameter to True if there are annotated disulfide bonds in the search database that you want to avoid searching, or when there is a known error in the disulfide bond annotation in the database. Errors in UniProt where disulfide bonds are placed on non-cysteine residues have been observed and can cause the search to crash.
Maximum Mass to Include PTMs	Limits the size of proteoforms that are considered with PTMs. It is unlikely to observe very large proteoforms, therefore it is unnecessary to consider all of the modified forms beyond the size expected to be observed in the experiment. Setting the limit to a reasonable size (100 kDa) allows the search to complete sooner.
N-Term Modifications to Include	N-terminal modifications are added at search time unless specifically added to the database. These modifications generally include acetylation or formulation, and the default is set to include N-terminal acetylation.
Delta M Mode	To identify unknown mass shifts, consider setting Delta M Mode to False. Searching with Delta M Mode will slow the search. Thermo Fisher Scientific recommends turning off FDR when utilizing Delta M Mode as the Delta M search affects the FDR threshold. Delta M Mode is suggested for data that has already been searched or in a target case.

Table 2. Advanced parameters for the Annotated Proteoform

Parameter	Description
Decoy Reps	Sets the number of times the decoy database is searched. The default setting is 1, which results in the fastest search, however results can vary minimally run-to-run: 1-2 proteoform differences base on exactly where the FDR thresholds are set. Increasing the number of times the decoy database is searched to 3 results in slower searches but reduces the amount of variation between analyses.
Maximum PrSMs Per Precursor	Limits the number of spectral matches attributed to a single precursor. Any given precursor can have many associated spectral matches depending on search parameters and the data. In general, there are very few highly confident matches (based on P-Score or E-Value.) This parameter dictates the top N scoring PrSMs to return per precursor. The default value is 3, however not all 3 PrSMs can pass the FDR cutoff depending on their confidence. Due to the nature of the MS isolation window, the occurrence of off-by-one Dalton deconvolution errors and modifications such as deamination, which are ~1 Dalton, do not expect to see multiple high-scoring PrSMs for a single precursor (this is why the default setting is 3). If you expect multiple proteoforms to have the same precursor mass, (for example, positional isomers such as histones), you can increase this setting to 5 or 10.
Minimum Matched Fragments	Lets you discard spectra that deconvoluted with less than <i>n</i> fragment ions. The default value is 3 and you can increase it to 5 or 10 without major loss of IDs. For smaller proteoforms (< 5kDa), set the minimum to 5 or lower. By discarding sparse spectra, the search completes faster.
Include Labile Modification Mass Shift	Lets you apply a mass shift to all precursors to account for loss of a labile modification. To activate the setting, set Include Labile Modification Mass Shift to True and input a mass in the Labile Modification Mass Shift parameter.

Perform a Subsequence search in the processing workflow

The Subsequence search comprehensively searches subsequences of the isoforms in a specified database. The node considers any PTMs or other annotations found in the database. The Subsequence search is more time- and resource-intensive than the Annotated Proteoform search. For any given sequence in the database, there could be hundreds of possible subsequences with additional complexity when considering PTMs and other modifications of each subsequence. Due to the increased search space, the Subsequence search utilizes more resources. Use the Subsequence search node to identify unexpected proteoforms in a sample and for discovering protein clips and endogenous truncations.

Thermo Fisher Scientific recommends first running the Subsequence search in the most restrictive manner possible and then incrementally increasing the search space as needed. A wide tolerance Subsequence search can take several days to complete and can draw an unsustainable amount of resource. To restrict the Subsequence search space, adjust the following parameters:

- For Precursor Mass Tolerance, set this to 15 ppm or fewer for complex databases.
- For small or targeted searches, tolerances of 2.2 Da can be useful to account for off-by-one errors.

❖ To add a Subsequence search in the processing workflow

1. Open a workflow template that includes the ProSightPD Subsequence Search node.
 - a. If the workflow does not include the Subsequence Search node, drag and drop the Subsequence Search node into the workflow.
 - b. Connect the node by dragging the arrows from the upstream node to the Subsequence Search node.
2. Select the Subsequence Search node to view the parameters.
3. Set the Input Database by selecting from the dropdown menu.

Tip To create a new database, launch the Database Manager. To select a newly created database, you must close and reopen any open studies.

4. Set the Precursor Mass Tolerance.

Thermo Fisher Scientific recommends 10 ppm.

Note Due to the large search space covered by the Subsequence Search, it is strongly suggested setting the precursor mass tolerance no greater than 25 ppm. For smaller data files or databases, you can use wider tolerances without causing the search time to increase.

5. Set the Fragment Mass Tolerance.

Thermo Fisher Scientific recommends 10 ppm. Use the narrowest tolerance possible for your data. As the tolerance widens, the search time and false positive rate increases.

6. Set the FDR calculation.

The default is True, which means that the FDRs are calculated based on the threshold set in the FDR node in the consensus workflow. The FDR node in the consensus workflow is required. If your sample is not complex or you want to see all PrSMs, set this value to False to not calculate the FDR.

Note Results are automatically displayed at the specified FDR cutoff for each result level (PrSM, Proteoform, Isoform, and so on).

7. (Optional) Select a static modification to apply to a specific residue.

You can apply one only static modification in the search node. To set additional modification, use the Protein Annotator tool in the Database Manager.

Advanced parameters of the Subsequence search

The default parameters of the Subsequence search are designed to work in most situations, however, most analyses can be improved/customized by optimizing the advanced parameters listed below. For a full list of advanced parameters, see “[Subsequence Search node.](#)” Table 3 is a list of advanced parameters and a brief description of how they impact the Subsequence search and results

Table 3. Advanced parameters of the Subsequence search

Parameter	Description
Maximum PTMs per Isoform	Limits the number of PTMs that are expanded and searched during search time. For example, if an isoform has 10 PTMs in the database and the maximum is set to 4, then only proteoforms with up to 4 PTMs are searched. This improves the search time. However, in cases of highly modified proteoforms, it is advised that you increase this value.
Maximum SNPs per Isoform	Limits the number of amino acid mutations of single nucleotide polymorphisms (SNPs) in the same way as in the Maximum PTMs per Isoform parameter.
Maximum PTMs per Proteoform	Limits the number of PTMs that are allowed on each proteoform during search time. For example, if an isoform has 10 PTMs in the database, the maximum per isoform is set to 8, and the maximum per proteoform is set to 4, then only proteoforms (or subsequences) with up to 4 PTMs are searched. This improves the search time. However, in cases of highly modified proteoforms, it is advised that you increase this value.

Tip The PSDB database contains isoforms and PTMs. All proteoforms to be searched are generated at search time according to the database and parameters selected in the search nodes. For the Subsequence search, each proteoform is expanded into its subsequences, so in the Subsequence search it is useful to limit the PTMs at the proteoform level.

Table 3. Advanced parameters of the Subsequence search

Parameter	Description
Maximum SNPs per Proteoform	Limits the number of amino acid mutations of single nucleotide polymorphisms (SNPs) in the same way as the Maximum PTMs per Proteoform parameter.
Remove Disulfide Bonds	Thermo Fisher Scientific recommends setting this parameter to False. You can set this parameter to True if there are annotated disulfide bonds in the search database that you want to avoid searching or when there is a known error in the disulfide bond annotation in the database. Errors in UniProt where disulfide bonds are placed on non-cysteine residues have been observed and can cause the search to crash.

Perform Single Proteoform searches in the processing workflow

The ProSightPD application uses the RESID™ database. The ProSightPD Single Proteoform Search node references all post-translational modifications as a truncated form of their RESID identification number.

❖ To perform single proteoform searches in the processing workflow

1. Open a workflow template that includes the ProSightPD Single Proteoform Search node, or add the node after the cRAWler node.
2. In a web browser, navigate to <https://proteininformationresource.org/cgi-bin/resid.pl>
3. Look up the RESID modification ID, and copy the integer that follows the “AA00.”

In this example, the integer is 76.

RESID Database Release 76.00 31-May-2018 List All Entries List All SequenceSpecs
Search version Beta-02

Sub-string text search for RESID ID, Name, Author, Citation, Title, Keyword, Feature and Enzyme; Exact match for SequenceSpec, CAS_ID, GO_ID, PDB_ID, PMID

Search (RESID ID ▼ 0076) and (RESID ID ▼) and Weight F chem ▼ ± 1.0 Submit


1 record found. (Note: In the following table, Fc=Weight F. chem.; Fp=Weight F. phys.; Cc=Weight C. chem.; Cp=Weight C. phys.)

RESID ID	Name	SequenceSpec	Weight	Keyword	Feature
AA0076	N6-methyl-L-lysine	L-lysine	Fc=142.20, Fp=142.110613, Cc=14.03, Cp=14.015650,	methylated amino acid	SP:MOD_RES N6-methylated lysine SP:MOD_RES N6-methyllysine

4. In the application, click the **ProSightPD Single Proteoform Search** node, and then click **Advanced Parameters**.
5. In the RESID Sequence, enter the protein sequence (Figure 17).

The RESID Sequence is the protein sequence plus the modification in the form of a RESID database entry.

Figure 17. The RESID Sequence entry for methylated lysine



RESID number in parentheses

Note Put the identifier in parenthesis before the amino acid containing the modification. For example, AA0049 N-acetyl-L-methionine becomes 49, and an acetylation of a methionine residue in a protein sequence is expressed as (49)M.

6. To add the modification, enter the RESID number in parentheses before the residue.
7. For manual precursor mass assignment, in the crawler node's Advanced Parameters, do the following:
 - a. Set Use Manual Precursor to **Always** or **If Not Detected**.
 - b. Type the expected Precursor Mass (Daltons).
8. Configure the consensus workflow, or click **Run**.

For basic workflow information, refer to the *Proteome Discoverer User Guide*.

For more information about using the node, see [Single Proteoform Search node](#).

Working with results

Contents


- The results report
- Quantitative results in ProSightPD
- Visualizing top-down results
- Filter results
- Visualize results with charts and graphs

The results report

The application generates a results report that contains a list of the matching proteoforms identified by the search nodes that you selected ([Figure 18](#)). You can access the results report through the application's Start page or File menu. Through the Job Queue, you can access the results, including a failed workflow analysis. The application can save multiple search results within a single report.

The application saves the results report in .pdResult file format so you can also access the results through the Windows Explorer.

Figure 18. Results reports displayed in the application window



The results report includes the following pages:

- [Proteins page](#)
- [Isoforms page](#)
- [PrSMs page](#)
- [MS/MS Spectrum Info page](#)
- [Input Files page](#)
- [Specialized Traces page](#)

Proteins page

By default the results report (.pdResult file) opens to the Protein Groups page.

Figure 19. Commonly used columns on the Proteins page

Proteins										
	Checked	Accession	# of Isoforms	# of Isoforms with Characterized Proteoforms	# of Proteoforms	Proteoform Ch	# Characterized Proteoforms	# of PrSMs	Q-value	
1	<input type="checkbox"/>	P64463	1	1	2	●	1	9	5.2E-039	
2	<input type="checkbox"/>	P33219	1	1	1	●	1	18	2.7E-035	
3	<input type="checkbox"/>	P0AB14	1	1	1	●	1	2	6.1E-033	
4	<input type="checkbox"/>	P0ACE7	1	1	1	●	1	6	8.7E-027	
5	<input type="checkbox"/>	P0A800	1	1	1	●	1	6	4.7E-025	
6	<input type="checkbox"/>	P0AFH8	6	5	9	●	7	98	2.2E-024	
7	<input type="checkbox"/>	P0AE08	1	1	6	●	2	26	2.4E-022	
8	<input type="checkbox"/>	P0A763	2	2	2	●	2	6	3.7E-021	
9	<input type="checkbox"/>	P0AC62	1	1	1	●	1	8	1.0E-020	
10	<input type="checkbox"/>	P0ADU5	1	1	1	●	1	3	1.2E-018	
11	<input type="checkbox"/>	P0AF50	3	2	4	●	2	5	1.9E-009	
12	<input type="checkbox"/>	P45470	1	1	1	●	1	2	1.2E-008	
13	<input type="checkbox"/>	P76402	1	1	1	●	1	2	1.3E-007	
14	<input type="checkbox"/>	P0ADE6	1	1	1	●	1	1	1.5E-007	
15	<input type="checkbox"/>	P45578	1	1	1	●	1	1	1.1E-004	

The Proteins page reports all of the proteins identified in the experiment that pass any user-defined threshold such as FDR. From this page, you can see the UniProt Accession and the number of proteoforms and isoforms of each protein identified in the experiment. Additionally, the number of PrSMs and characterized proteoforms (that is, proteoforms with C-Score greater than 40) are listed. If you enabled FDR, then the Q-Value for each protein is also displayed.

You can access additional information for each protein identified by selecting it and clicking the down arrow next to Show Associated Tables ().

Table 4. Proteins page columns

Column	Description
Accession	UniProt™ accession identifiers used to cite UniProt entries.
# of Isoforms	The number of isoforms identified from the associated protein accession.
# of Isoforms with Characterized Proteoforms	The number of isoforms with proteoforms that have a C-Score greater than 40.
# Proteoforms	The number of identified proteoforms from the associated UniProt protein accession
Proteoform Characterization Confidence	The visual indication of C-Score confidence (>40=green,<40->3=yellow, <3=red). You can define these thresholds in the ProSight Protein and Proteoform Filter node in the consensus workflow.
# Characterized Proteoforms	The number of proteoforms having a C-Score greater than 40.
Q-value	The likelihood of discovering the associated protein by chance given the experimental evidence.

Isoforms page

This page displays all of the identified isoforms (Figure 20). The page also includes the UniProt description for each isoform and the mass in kDa, pI, and, if enabled, the protein annotations for gene ontologies and various pathway analysis.

You can access additional information for each isoform identified by selecting the isoform and clicking the down arrow next to Show Associated Tables (Show Associated Tables).

Figure 20. The Isoforms page

	Checked	Proteoform Characterization Confidence	Accession	Description	# Proteoforms	# PrSMs	# AAs	MW [kDa]	calc. pI	Q-value	# Characterized Proteoforms
1	<input type="checkbox"/>	●	P60723	50S ribosomal protein L4	2	6	201	22.1	9.73	4.0E-004	2
2	<input type="checkbox"/>	●	P0AFX0	Ribosome hibernation promoting factor	3	3	95	10.7	7.05	9.3E-005	1
3	<input type="checkbox"/>	●	P0A7W1	30S ribosomal protein S5	2	2	167	17.6	10.11	2.1E-012	1
4	<input type="checkbox"/>	●	P0A7S9	30S ribosomal protein S13	2	27	118	13.1	10.78	2.3E-031	2
5	<input type="checkbox"/>	●	P02413	50S ribosomal protein L15	1	3	144	15.0	11.18	8.2E-008	1
6	<input type="checkbox"/>	●	P60438	50S ribosomal protein L3	3	8	209	22.2	9.91	6.3E-012	2
7	<input type="checkbox"/>	●	P0A7L8	50S ribosomal protein L27	1	14	85	9.1	10.58	7.1E-010	1
8	<input type="checkbox"/>	●	P0AG59	30S ribosomal protein S14	1	14	101	11.6	11.17	5.9E-016	1
9	<input type="checkbox"/>	●	P68919	50S ribosomal protein L25	3	11	94	10.7	9.60	7.4E-050	1
10	<input type="checkbox"/>	●	P0A7W7	30S ribosomal protein S8	1	10	130	14.1	9.42	4.6E-012	1
11	<input type="checkbox"/>	●	P0A7J3	50S ribosomal protein L10	3	48	165	17.7	8.98	4.1E-057	2
12	<input type="checkbox"/>	●	P0A7L3	50S ribosomal protein L20	1	28	118	13.5	11.47	3.1E-042	1
13	<input type="checkbox"/>	●	P0A7R9	30S ribosomal protein S11	1	13	129	13.8	11.33	2.7E-023	1
14	<input type="checkbox"/>	●	P02359	30S ribosomal protein S7	1	14	179	20.0	10.36	4.3E-057	1
15	<input type="checkbox"/>	●	P0A7U3	30S ribosomal protein S19	1	19	92	10.4	10.52	2.3E-078	1
16	<input type="checkbox"/>	●	P68679	30S ribosomal protein S21	1	11	71	8.5	11.15	3.2E-031	1
17	<input type="checkbox"/>	●	P0ADZ0	50S ribosomal protein L23	1	10	100	11.2	9.94	1.1E-061	1
18	<input type="checkbox"/>	●	P0A7T7	30S ribosomal protein S18	1	9	75	9.0	10.59	1.3E-045	1
19	<input type="checkbox"/>	●	P0C018	50S ribosomal protein L18	1	3	117	12.8	10.42	4.4E-017	1
20	<input type="checkbox"/>	●	P0ADZ4	30S ribosomal protein S15	1	8	89	10.3	10.40	1.3E-022	1
21	<input type="checkbox"/>	●	P0A7V8	30S ribosomal protein S4	1	5	206	23.5	10.05	3.3E-027	1
22	<input type="checkbox"/>	●	P0AG55	50S ribosomal protein L6	2	6	177	18.9	9.70	1.3E-027	1
23	<input type="checkbox"/>	●	P0A7U7	30S ribosomal protein S20	1	16	87	9.7	11.18	2.8E-022	1
24	<input type="checkbox"/>	●	P60624	50S ribosomal protein L24	3	27	104	11.3	10.21	1.2E-085	3

Table 5. Commonly used columns on the Isoforms page (Sheet 1 of 2)

Column	Description
Proteoform Characterization Confidence	Visual indication of C-Score confidence (>40=green,<40->3=yellow, <3=red). For information on defining the thresholds, see the PrSM Grouper node parameters.
Accession	UniProt accession identifiers used to cite UniProt entries.
Description	The protein description from UniProt.
# Proteoforms	The number of identified proteoforms from the associated UniProt protein accession.
# PrSMs	The number of PrSMs supporting this isoform assignment.

Table 5. Commonly used columns on the Isoforms page (Sheet 2 of 2)

Column	Description
# Characterized Proteoforms	The number of proteoforms having a C-Score greater than 40.
Q-value	The chance of discovering the associated result given the experimental evidence.

Proteoforms page

The Proteoforms page displays all of the identified proteoforms.

Figure 21. The Proteoforms page

The Proteoforms page also includes the identified sequence for each proteoform and the number of supporting PrSMs, theoretical mass in kDa, search engine confidence, -Log P-Score, -Log E-Value, best PrSM C-score, Sequence Coverage (as % residue cleavage.) If FDR is enabled, proteoform level Q-Values are also displayed.

For each result at the proteoform level, you can launch the ProSight Lite and TDValidator Lite applications. For additional information on these applications, see:

- [Use ProSight Lite](#)
- [Use TDValidator Lite](#)



You can access additional information for each identified proteoform by selecting the proteoform and clicking the down arrow next to Show Associated Tables ( Show Associated Tables).

Table 6. Commonly used columns on the Proteoforms page

Column	Description
Proteoform Characterization Confidence	Visual indication of C-Score confidence (>40=green,<40->3=yellow, <3=red). For information on defining the parameters, see PrSM Grouper node .
Sequence	Proteoform sequence.
# PrSMs	Number of proteoforms.
Theo Mass [Da]	Proteoform mass based on the matched sequence.
Confidence	Visual representation of the proteoform's PrSM confidence based on the -log E-value levels specified in the ProSightPD configuration. Higher confidence = green, medium confidence = yellow, and low confidence = red
-Log P-Score	Statistical confidence for proteoform assignment. Large positive numbers are more confident. Values less than 4 are considered poor.
-Log E-Value	Statistical confidence for proteoform assignment adjusted for the search space.
Best PrSM C-Score	Highest C-Score of PrSMs supporting the associated proteoform.
Average PrSM Detected Neutral Mass	Average observed mass from all PrSMs associated with a proteoform.
% Residue Cleavages	Percentage of observed fragments compared to the total possible cleavages. Also known as sequence coverage.
External Top Down Displays	Click to view the result in the ProSight Lite or the TDValidator Lite applications. For more information, see Use ProSight Lite and Use TDValidator Lite .

PrSMs page

The Proteoform Spectral Matches (PrSMs) page reports all of the PrSMs identified in the experiment.

The annotated sequence showing modifications, LogE,-LogP, are shown for each PrSM. To view additional information about related proteoforms, proteins, and isoforms, click the down arrow next to Show Associated Tables ( Show Associated Tables).

When Delta M Mode is enabled, fragments found using the Delta M Mode display on the PrSMs page, along with two additional columns: Delta M Residue and Delta M Fragment(s), where the delta mass and the possible residue where the delta mass was localized. The delta masses also display in the ProSightPD fragment map with an orange square around the potential residue where the delta mass was localized.

Figure 22. PrSMs page

Confidence	Protein Accessions	Annotated Sequence	Modifications	Theo. Mass [Da]	Mass [Da]	Δ Mass [ppm]	Δ Mass [Da]	Corrected Delta Mass (ppm)	Corrected Delta Mass (Da)	Q-value	Delta M Residue	Delta M Fragment(s)
1	P02358	MRHYEIFMVHPDC	I71(DeltaM 12)	15301.53173	15301.53150	-0.02	-0.00023	0.00	0.000	7.3E-006	I71 (124.0924753)	Y110
2	P02358	MRHYEIFMVHPDC		15177.43925	15291.48380	7514.08	114.04455	0.00	0.000	1.7E-005		
3	P02358	MRHYEIFMVHPDC		15177.43925	15317.49230	9227.71	140.05305	9.227.73	140.053	1.7E-005		
4	P02358	MRHYEIFMVHPDC	E125(DeltaM)	15163.46353	15163.46330	-0.02	-0.00023	0.00	0.000	1.4E-012	E125 (-13.975724)	Y31, Y49, Y55, Y56, Y
5	P02358	MRHYEIFMVHPDC	K93(N6-acety)	15292.55993	15292.55970	-0.02	-0.00023	0.00	0.000	2.8E-009	D114 (73.110110)	Y49, Y56, Y56, Y57, Y
6	P02358	MRHYEIFMVHPDC	D114(DeltaM)	15292.55993	15292.55970	-0.02	-0.00023	0.00	0.000	2.8E-009	D114 (115.12067)	Y49, Y56, Y56, Y57, Y
7	P02358	MRHYEIFMVHPDC		15177.43925	15301.50870	8174.60	124.06945	0.00	0.000	6.0E-006		
8	P02358	MRHYEIFMVHPDC		15177.43925	15316.53780	9164.82	139.09855	0.00	0.000	6.0E-006		
9	P02359	PRRRVIGQRKLPDI		17462.40726	17462.48610	4.51	0.07884	0.00	0.000	1.9E-014		
10	P02359	PRRRVIGQRKLPDI		17462.40726	17445.43170	-972.12	-16.97556	0.00	0.000	1.8E-014		
11	P02359	PRRRVIGQRKLPDI		17462.40726	17462.53010	7.03	0.12284	0.00	0.000	9.6E-015		
12	P02359	PRRRVIGQRKLPDI	L123(DeltaM)	17462.53042	17462.53010	-0.02	-0.00032	0.00	0.000	3.8E-013	L123 (0.12315531)	B123, B124, Y45, Y66
13	P02359	PRRRVIGQRKLPDI		17462.40726	17462.48610	4.51	0.07884	4.53	0.079	7.9E-018		
14	P02359	PRRRVIGQRKLPDI		17462.40726	17462.53010	7.03	0.12284	7.05	0.123	3.2E-018		
15	P02359	PRRRVIGQRKLPDI		17462.40726	17462.53010	7.03	0.12284	7.05	0.123	4.8E-016		
16	P02413	MRLNTLSPAEGSK	E86(DeltaM)	14957.33604	14957.33570	-0.02	-0.00034	0.00	0.000	3.0E-020	E86 (-13.9874546)	B89, B92, B95, B99, E
17	P02413	MRLNTLSPAEGSK	E86(DeltaM)	14957.38464	14957.38430	-0.02	-0.00034	0.00	0.000	3.3E-009	E86 (-13.9388546)	B86, B98, B100, B102
18	P02413	MRLNTLSPAEGSK	D90(DeltaM 1)	14961.33784	14961.33750	-0.02	-0.00034	0.00	0.000	6.4E-004	D90 (121.054835)	B135, Y67, Y135, Y13
19	P02413	MRLNTLSPAEGSK		14971.32349	14996.27320	1666.50	24.94971	1.666.52	24.950	1.7E-003		
20	P02413	MRLNTLSPAEGSK		14971.32349	14973.37480	137.02	2.05131	0.00	0.000	1.7E-003		
21	P02413	MRLNTLSPAEGSK		14971.32349	14941.41350	-1997.82	-29.90999	-1.997.80	-29.910	6.2E-003		
22	P02413	MRLNTLSPAEGSK		14971.32349	14900.28200	-4745.17	-71.04149	0.00	0.000	6.2E-003		
23	P02413	MRLNTLSPAEGSK	R48(DeltaM -1)	14912.29724	14912.29690	-0.02	-0.00034	0.00	0.000	6.2E-003	R48 (-59.0262546)	B136
24	P0A6K3	AKGQSLQDFLNAL		11028.54188	11034.58772	548.20	6.04584	0.00	0.000	2.2E-008		
25	P0A707	mKGGKRVQTRPN	M1(N-methyl-)	20565.06894	20564.11880	-46.20	-0.95014	0.00	0.000	8.8E-006		
26	P0A7J3	ALNLQDKQAIAVEV	K104(N6-acet)	17569.38507	17569.38470	-0.02	-0.00037	0.00	0.000	2.5E-011	E113 (-41.975509)	B113, B115, B116, B1
27	P0A7J3	ALNLQDKQAIAVEV	K36(N6-acety)	17569.38507	17569.38470	-0.02	-0.00037	0.00	0.000	2.5E-011	R60 (-41.9755096)	B75, B77, B79, B80, E
28	P0A7J3	ALNLQDKQAIAVEV		17569.35002	17562.44670	-962.09	-16.90332	0.00	0.000	3.9E-010		
29	P0A7J3	ALNLQDKQAIAVEV	K36(N6-acety)	17568.32067	17568.32030	-0.02	-0.00037	0.00	0.000	6.2E-003	E115 (-85.050474)	B115, Y105
30	P0A7J3	ALNLQDKQAIAVEV	K36(N6-acety)	17568.32067	17568.32030	-0.02	-0.00037	0.00	0.000	6.2E-003	E115 (-43.039909)	B115, Y105
31	P0A7J3	ALNLQDKQAIAVEV		17569.35002	17568.34760	-57.05	-1.02424	0.00	0.000	6.2E-003		
32	P0A7J3	ALNLQDKQAIAVEV		17569.35002	17569.38470	1.97	0.03468	2.00	0.035	5.1E-013		
33	P0A7J3	PASQIDRLATLPITYE		4997.67943	4997.68250	0.61	0.00307	0.63	0.003	2.9E-029		
34	P0A7J3	ALNLQDKQAIAVEV		16216.59757	16216.62400	1.63	0.02643	1.65	0.027	1.3E-010		

To further investigate each result, you can launch the ProSight Lite and the TDValidator Lite applications from the PrSMs table.

Table 7. Commonly used columns on the PrSMs page (Sheet 1 of 2)

Column	Description
Confidence	Visual representation of confidence. High Confidence = green, Medium=Yellow, Low=Red. You can define these thresholds.
Detected Ion Count	The number of fragment ions found in the fragment ion spectra.
Identifying Node	Search node used to identify this PrSM.
Annotated Sequence	The proteoforms backbone sequence and PTMs.
Modifications	List of modifications found on each proteoform.
# Proteins	The number of proteins each proteoform maps back to.
Protein Accessions	Parent protein accession.
Charge	Charge state isolated for fragmentation of the proteoform.
m/z [Da]	Observed mass.
Mass [Da]	Monoisotopic observed Mass.
Theo. Mass [Da]	Theoretical mass based on matched proteoform from the database.
Δ Mass [ppm]	Observed theoretical mass error in parts per million.

Table 7. Commonly used columns on the PrSMs page (Sheet 2 of 2)

Column	Description
Δ Mass [Da]	Observed-theoretical mass error in Daltons.
Matched Ions	Number of matched fragment ions supporting the PrSM.
Delta M Residue	Possible location of the delta mass supported by the corresponding Delta M fragments.
Delta M Fragment(s)	Fragments found using Delta M Mode.


MS/MS Spectrum Info page

The MS/MS Spectrum Info page lists all of the MS/MS spectra searched in the experiment and certain details such as scan number, injection time, # of PrSM (supported by this MS/MS scan), # of Precursors (associated with this MS/MS scan). You can use the MS/MS Spectrum Info table to interrogate which MS/MS spectra are associated with which proteoform or PrSM result.

Figure 23. MS/MS Spectrum Info page

Proteins		Isoforms		Proteoforms		PrSMs		MS/MS Spectrum Info		Input Files		Specialized Traces		Study Information	
	File	RT [min]	First Scan	Mass Analyzer	Activation Type	NCE [%]	MS Order	# PrSMs	Isolation Interference [%]	Ion Inject Time [ms]	# Precursors	# Identified Precursors	Precursor m/z [Da]	Precursor MH+ [Da]	Precursor C
1	F2	47.0211	2977	FTMS	HCD	20.0	MS2	0	92	3.652	1	0	100.40984	895.63038	
2	F2	47.0211	2977	FTMS	HCD	20.0	MS2	1	92	3.652	1	0	2161.52152	19445.63548	
3	F2	47.0211	2977	FTMS	HCD	20.0	MS2	0	92	3.652	1	0	113.08955	1009.74778	
4	F2	47.0211	2977	FTMS	HCD	20.0	MS2	0	92	3.652	1	0	933.02067	8389.12778	
5	F2	47.0211	2977	FTMS	HCD	20.0	MS2	0	92	3.652	1	0	865.26018	18150.31818	
6	F2	47.0254	2978	FTMS	HCD	20.0	MS2	3	68	7.616	1	0	1009.30233	18150.31818	
7	F2	47.0254	2978	FTMS	HCD	20.0	MS2	0	68	7.616	1	0	1009.74778	1009.74778	
8	F2	47.0254	2978	FTMS	HCD	20.0	MS2	0	68	7.616	1	0	19445.63548	19445.63548	
9	F2	47.0254	2978	FTMS	HCD	20.0	MS2	0	68	7.616	1	0	895.63038	895.63038	
10	F2	47.0254	2978	FTMS	HCD	20.0	MS2	1	68	7.616	1	0	467.01397	8389.12778	
11	F2	47.0298	2979	FTMS	HCD	20.0	MS2	3	86	3.381	1	0	908.47282	18150.31818	
12	F2	47.0298	2979	FTMS	HCD	20.0	MS2	1	86	3.381	1	0	420.41330	8389.12778	
13	F2	47.0298	2979	FTMS	HCD	20.0	MS2	0	86	3.381	1	0	51.44430	1009.74778	
14	F2	47.0298	2979	FTMS	HCD	20.0	MS2	0	86	3.381	1	0	973.23869	19445.63548	
15	F2	47.0298	2979	FTMS	HCD	20.0	MS2	0	86	3.381	1	0	45.73843	895.63038	
16	F2	47.0342	2980	FTMS	HCD	20.0	MS2	3	87	5.669	1	0	865.26018	18150.31818	
17	F2	47.0342	2980	FTMS	HCD	20.0	MS2	1	87	5.669	1	0	400.44159	8389.12778	
18	F2	47.0342	2980	FTMS	HCD	20.0	MS2	0	87	5.669	1	0	926.94195	19445.63548	
19	F2	47.0342	2980	FTMS	HCD	20.0	MS2	0	87	5.669	1	0	43.60838	895.63038	
20	F2	47.0342	2980	FTMS	HCD	20.0	MS2	0	87	5.669	1	0	49.04254	1009.74778	
21	F2	47.0653	2982	FTMS	HCD	20.0	MS2	3	87	4.469	1	0	825.97519	18150.30148	
22	F2	47.0653	2982	FTMS	HCD	20.0	MS2	1	87	4.469	1	0	382.28562	8389.13078	
23	F2	47.0653	2982	FTMS	HCD	20.0	MS2	3	87	4.469	1	0	943.29919	20731.42948	
24	F2	47.0696	2983	FTMS	HCD	20.0	MS2	1	62	23.302	1	0	560.21551	8389.13078	
25	F2	47.0696	2983	FTMS	HCD	20.0	MS2	3	62	23.302	1	0	1383.03542	20731.42948	
26	F2	47.0696	2983	FTMS	HCD	20.0	MS2	2	62	23.302	1	0	1210.96022	18150.30148	
27	F2	47.0740	2984	FTMS	HCD	20.0	MS2	3	70	2.079	1	0	956.23329	18150.30148	
28	F2	47.0740	2984	FTMS	HCD	20.0	MS2	1	70	2.079	1	0	442.48746	8389.13078	
29	F2	47.0740	2984	FTMS	HCD	20.0	MS2	3	70	2.079	1	0	1092.08213	20731.42948	
30	F2	47.0784	2985	FTMS	HCD	20.0	MS2	3	59	4.071	1	0	1068.61282	18150.30148	
31	F2	47.0784	2985	FTMS	HCD	20.0	MS2	1	59	4.071	1	0	494.42631	8389.13078	
32	F2	47.0784	2985	FTMS	HCD	20.0	MS2	3	59	4.071	1	0	1220.44388	20731.42948	

❖ Work with the MS/MS Spectrum Info page

- To display the associated precursor scans (precursor before and after the MS/MS scan), click the **MS/MS Spectrum Info** icon () in the upper-left corner of the Field Chooser.

The Field Chooser menu opens.

- From the Field Chooser menu, select **Master Scan**.

- To display Associated PrSMs, Proteoforms, and so on, select the down arrow next to Show Associated Tables ( Show Associated Tables).

Table 8. Commonly used columns on the MS/MS Spectrum Info page

Column	Description
File	The .raw data file name.
RT [min]	Retention time of the MS/MS scan (in minutes).
First Scan	First precursor scan associated with this MS/MS scan.
Activation Type	Activation type used to generate ions in the MS/MS scan.
NCE[%]	Normalized Collision Energy applied if HCD is used.
# PrSMs	Number of PrSMs that this MS/MS scan supports.
Isolation Interference [%]	Measure of co-isolated precursor interference.
Ion Inject Time [ms]	How long ions were collected prior to analysis. Shorter times indicate a more abundant signal.
m/z [Da]	Observed mass.
Mass [Da]	Monoisotopic observed mass.
Theo. Mass [Da]	Theoretical mass based on matched proteoform from database.
Δ Mass [ppm]	Observed theoretical mass error in parts per million.
Δ Mass [Da]	Observed theoretical mass error in Daltons.
Matched Ions	Number of matched fragment ions supporting the PrSM.


Input Files page

The Input Files page shows the raw data files searched in this experiment. Refer to the *Proteome Discoverer User Guide* for more information.

Specialized Traces page

The Specialized Traces page provides the total ion chromatogram (TIC) and base peak chromatogram (BPC) for each input file ([Figure 24](#)).

Figure 24. Specialized Traces page



	Checked	Tags	File ID	Trace/Detector Type	Description	Scan Filter	Spectrum File
1	<input type="checkbox"/>	○○○○○	F1	TIC Trace	MS Order: MS1	FTMS + p NSI Full ms [600.0000-2000.0000]	SEC_Fraction 14.raw
2	<input type="checkbox"/>	○○○○○	F1	BPC Trace	MS Order: MS1	FTMS + p NSI Full ms [600.0000-2000.0000]	SEC_Fraction 14.raw

[Table 9](#) lists the columns of the Specialized Traces page.

Table 9. Columns of the Specialized Traces page

Column	Description
Tags	Not used with the ProSightPD workflows supplied with the software.
File ID	Number the application assigns when importing the raw data file.
Trace/Detector Type	Base Peak Chromatogram (BPC) or Total Ion Chromatogram (TIC).
Description	MS or MS/MS.
Scan Filter	The <i>m/z</i> range.
Spectrum File	The raw data file name.

Quantitative results in ProSightPD

The LFQ workflow generates additional results pages and columns:

- [Feature Group page](#)
- [Consensus Feature Group page](#)

Feature Group ID page

The Feature Group page reports a list of all features individually measured from all data files. Use the Feature Group page to see the mass and abundance of each detected feature as well as additional features (Figure 25).

Figure 25. Feature Group page columns

	Checked	File ID	Apex RT [min]	Area	Intensity	# LCMS Features	Scan Range	RT Range	Mass [Da]	Quan Info	# PrSMs
1	<input type="checkbox"/>	F44	61.74	5.108e6	0.000e0	1	1595-1614	61.383-62.185	1996.14556	Shared	2
2	<input type="checkbox"/>	F44	19.00	1.620e6	0.000e0	1	505-515	18.828-19.203	2011.19391	Shared	2
3	<input type="checkbox"/>	F44	33.48	8.849e7	0.000e0	1	889-919	32.999-33.91	2012.22612	Shared	4
4	<input type="checkbox"/>	F44	89.62	4.979e6	0.000e0	1	2231-2249	89.418-90.253	2014.12513	Shared	3
5	<input type="checkbox"/>	F44	16.99	2.302e6	0.000e0	1	445-458	16.696-17.181	2029.14486		0
6	<input type="checkbox"/>	F44	17.69	5.041e7	0.000e0	1	459-494	17.22-18.369	2038.08733	Shared	6
7	<input type="checkbox"/>	F44	39.52	1.788e6	0.000e0	1	1060-1074	39.253-39.768	2047.00339	Shared	2
8	<input type="checkbox"/>	F44	37.82	1.987e7	0.000e0	1	1010-1029	37.46-38.14	2055.21861	Shared	2
9	<input type="checkbox"/>	F44	25.36	3.569e6	0.000e0	1	670-679	25.209-25.535	2057.17140	Shared	2
10	<input type="checkbox"/>	F44	20.35	5.146e7	0.000e0	1	535-559	19.904-20.881	2070.18950		1
11	<input type="checkbox"/>	F44	25.44	4.166e7	0.000e0	1	665-690	25.003-25.911	2075.12476	Shared	2
12	<input type="checkbox"/>	F44	55.10	2.843e5	0.000e0	1	1440-1449	54.89-55.281	2077.20110		0
13	<input type="checkbox"/>	F44	35.34	8.595e7	0.000e0	1	945-974	34.889-36.153	2083.26645	Shared	2
14	<input type="checkbox"/>	F44	95.62	2.173e6	0.000e0	1	2361-2379	95.412-96.25	2085.16189	Shared	2
15	<input type="checkbox"/>	F44	15.40	9.111e6	0.000e0	1	405-419	15.137-15.713	2091.99567		0
16	<input type="checkbox"/>	F44	22.87	1.410e6	0.000e0	1	605-614	22.7-23.042	2099.20770		0
17	<input type="checkbox"/>	F44	23.27	5.887e6	0.000e0	1	610-629	22.885-23.648	2105.10250		2

Table 10 lists the Feature Group page columns.

Table 10. Feature Group page columns

Column	Description
File ID	A short identifier of the input file in which the application detected the feature.
Apex RT [min]	The retention time (RT) at which the abundance is greatest for the feature.
Area	The measured peak area for a feature.
Intensity	Displays the peak height of the largest theoretical contained LCMS peak.
# LCMS Features	The number of LC/MS features found within your defined mass and retention time thresholds. This value is usually one. When many feature groups contain more than 2 or 3 LC/MS features, consider tightening the RT and Mass tolerance in the feature detector node. For information on feature detector nodes and their parameters, see Feature Detector nodes .
Scan Range	The scan numbers where the feature was detected.
RT Range	The retention time range when the feature was detected.
Mass (Da)	The detected mass of the feature.
Quan Info	The column indicates if the feature is shared or unique. <ul style="list-style-type: none"> • <i>Shared</i> indicates features mapping to multiple PrSMs • Empty column (indicating unique) indicates that the feature maps to 1 or 0 PrSMs
#PrSMs	The number of PrSMs.

Consensus Feature Group page

Consensus feature groups are a collection of similar (within mass and RT tolerance) feature groups, one from each file, representative of the feature across files. The Consensus Feature Groups page provides a list of all the consensus feature groups found.

[Figure 26](#) shows some of the columns of the Consensus Feature Group page.

Figure 26. Consensus Feature Group page (detail)

		Proteins	Isoforms	Proteoforms	PrSMs	MS/MS Spectrum Info	Input Files	Specialized Traces	Consensus Feature Group	Feature Group					
	Checked	Avg. Apex RT [min]	Mass [Da]	# Files	Abundances (Normalized)										Abundance
1	<input type="checkbox"/>	0.53	5427.92083	9	9.806e6		9.590e6	2.465e7	1.564e7	3.228e6	1.474e7	1.015e7	6.825e6	4.380e6	3.938e6
2	<input type="checkbox"/>	0.17	9376.18216	5						2.205e6	1.506e7	1.853e7	1.070e7	8.275e6	
3	<input type="checkbox"/>	0.47	7028.43143	5						2.836e7	4.402e7	4.918e7	2.223e7	2.182e7	
4	<input type="checkbox"/>	0.30	9234.27492	5						8.096e5	6.060e6	7.454e5	6.855e5	5.529e5	
5	<input type="checkbox"/>	0.22	4090.23993	1						6.638e6					
6	<input type="checkbox"/>	0.67	10363.19601	5						2.361e7	1.404e7	1.789e7	9.918e6	4.614e6	
7	<input type="checkbox"/>	0.35	9234.07360	5						7.104e5	6.060e6	7.454e5	6.855e5	5.529e5	
8	<input type="checkbox"/>	0.64	10348.58096	5						6.280e6	1.682e7	1.216e7	7.124e6	1.216e7	
9	<input type="checkbox"/>	0.97	5403.15377	4						2.377e6		8.987e5	1.127e6	1.140e6	
10	<input type="checkbox"/>	1.40	5954.98406	10	1.798e7	2.030e7	1.337e7	2.514e7	2.593e7	1.060e8	1.237e8	7.262e7	8.975e7	5.969e7	7.222e6
11	<input type="checkbox"/>	1.30	6007.89833	4						2.106e6	2.289e6		1.571e6	1.787e6	
12	<input type="checkbox"/>	1.33	9634.23651	5						1.438e6	3.380e6	2.950e6	2.155e6	1.600e6	
13	<input type="checkbox"/>	1.26	9647.25005	5						1.567e6	3.637e6	1.625e6	6.141e5	2.915e6	
14	<input type="checkbox"/>	1.43	10094.74572	10	6.071e6	6.677e6	6.102e6	1.726e7	1.178e7	1.728e6	2.902e6	2.605e6	1.765e6	1.220e6	2.438e6
15	<input type="checkbox"/>	1.45	10175.21336	10	8.058e6	1.611e7	4.231e6	2.061e7	1.794e7	7.272e6	2.334e6	6.284e6	6.158e6	3.097e6	3.236e6

Table 11 lists the columns on the Consensus Feature Group page.

Table 11. Consensus Feature Group page columns

Column	Description
Average Apex RT	The average apex retention time of all of the feature groups composing the consensus feature group.
Mass [Da]	The average mass of the feature groups composing the consensus feature groups.
# Files	The number of files where a feature was detected.
Abundance (Normalized)	The normalized abundance of each feature group in the consensus feature group.
Abundances (per File)	The Abundances column reports the measured abundance of each feature group in the consensus feature group.
Quan Info	<i>Shared</i> indicates that a consensus feature maps back to two or more proteoforms. This occurs when there are two species of similar masses with similar RT, which can occur when there are two isobaric proteoforms of the same isoform.

Quantitative results columns in the Proteoforms page

When you use a label-free quantitation workflow, quantitative results appear in the Proteoforms Table.

For information on the other columns in the Proteoform page, see [Proteoforms page](#).

Figure 27. Quantitative information columns in the Proteoforms page

ms	Proteoforms	PrSMs	MS/MS Spectrum Info	Input Files	Specialized Traces	Consensus Feature Group	Feature Group						
Abundance Ratio	Abundance Ratios (log2)	Abundance Ratio Adj. P-Value	Abundances (Group)	Abundances									
69.098	6.11	5.5e-4	957290458.0	13854114.5	2.576e8	2.620e8	4.754e8	2.832e8	4.167e8	1.667e6	4.033e6	2.909e6	3.850e6
100.000	6.64	1.8e-3	039185949.4	6431996.2	4.173e8	3.542e8	7.657e8	3.112e8	4.389e8	2.860e6	7.728e5	2.163e6	
0.060	-4.06	1.8e-3	3572017.0	59685871.4			9.679e5	3.607e6		1.165e7	1.484e7	8.507e6	
48.250	5.59	1.8e-3	39156640.1	811536.8	9.466e6	1.169e7	2.741e7	8.063e6	1.700e7		2.079e5		
48.250	5.59	1.8e-3	39156640.1	811536.8	9.466e6	1.169e7	2.741e7	8.063e6	1.700e7		2.079e5		
0.100	-3.32	1.8e-3	569378.5	5677553.2	2.398e5	1.486e5	3.959e5		3.003e5		1.402e6	9.766e5	3.361e6
11.711	3.55	2.3e-3	106773027.2	9117141.3	3.897e7	2.671e7	1.168e8		7.244e7		2.888e6	2.071e6	3.749e6
35.491	5.15	2.4e-3	125522829.2	3536773.5	3.722e7	5.004e7	5.435e7	4.812e7	5.449e7	8.108e5	9.217e5	5.912e5	9.870e5
0.103	-3.27	2.4e-3	1529018.7	14790254.5	3.382e5	4.184e5	6.636e5	2.689e5	6.962e5		5.580e6	2.555e6	5.374e6
8.412	3.07	2.5e-3	56866413.9	6760098.8	2.049e7	1.556e7	3.458e7	1.737e7	2.792e7	1.515e6	1.718e6	8.194e5	3.017e6
0.045	-4.48	2.5e-3	173312.6	3868365.8			1.783e5		7.311e4	1.189e6	7.725e5	6.467e5	1.447e6
8.412	3.07	2.5e-3	56866413.9	6760098.8	2.049e7	1.556e7	3.458e7	1.737e7	2.792e7	1.515e6	1.718e6	8.194e5	3.017e6
0.226	-2.14	2.6e-3	20302705.5	89750018.7	7.222e6	5.556e6	1.337e7	4.083e6	1.126e7	1.936e7	3.143e7	1.214e7	4.893e7
35.727	5.16	2.6e-3	342332892.3	9581833.0	1.375e8	1.139e8	1.798e8	9.134e7	1.240e8	1.723e6			
28.686	4.84	2.6e-3	177073795.0	6172778.4	5.549e7	5.256e7	1.389e8	4.509e7	7.687e7	1.748e6	1.140e6	1.032e6	9.885e5
28.686	4.84	2.6e-3	177073795.0	6172778.4	5.549e7	5.256e7	1.389e8	4.509e7	7.687e7	1.748e6	1.140e6	1.032e6	9.885e5
35.727	5.16	2.6e-3	342332892.3	9581833.0	1.375e8	1.139e8	1.798e8	9.134e7	1.240e8	1.723e6			
0.179	-2.48	2.8e-3	786820.4	4400007.0	3.090e5	2.536e5	5.151e5		3.493e5	1.150e6	7.149e5	3.219e6	

Table 12 lists the columns that contain quantitative information. For the other columns, see Table 6.

Table 12. Quantitative information columns in the Proteoforms page


Column	Description
Abundance Ratio	The ratio of the grouped abundances based on requested ratios from the grouping and analysis tab.
Abundance Ratios (log2)	The abundance ratio as a logarithm.
Abundance Ratio Adjusted P-Value	The p-values adjusted by using the Benjamini-Hochberg correction for the false discovery rate. For more information on p-values, refer to Calculating P-Values and Adjusted P-Values for Quantification Results in the <i>Proteome Discoverer User Guide</i> .
Abundances Grouped	Displays the abundance values of the sample groups. A grouped abundance value is calculated as the arithmetic mean of all the replicate abundance values within a sample group. This column appears when you group samples in the analysis setup and there is at least one sample group consisting of at least two samples. You can specify the sample grouping on the Grouping and Quantification page when you set up an analysis.
Abundances	The abundance values.

Visualizing top-down results

These topics describe tools for viewing top-down results:

- [Use ProSight Lite](#)
- [Use TDViewer](#)
- [Use TDValidator Lite](#)


Use ProSight Lite

To display the sequence maps for selected proteoforms in the results view, click the **Show ProSightPD Fragment Map** icon, . These sequence maps are annotated with modifications and fragments from experimental MS/MS data.

If the ProSight Lite application is not already installed, install it from the prosightlite.northwestern.edu website.

By clicking PSLite in the External Top Down Displays column or from the ProSightPD Fragment Map, you can view the fragment map associated with the match in the ProSight Lite application ([Figure 28](#)).

Figure 28. Associated fragment map



Use TDViewer

ProSightPD nodes are integrated with the TDViewer application. The TDViewer application lets you share Proteome Discoverer results outside of the application.

❖ To install the TDViewer

1. To download the TDViewer, navigate to <http://topdownviewer.northwestern.edu/>.
2. Follow the viewer's download and installation instructions.

❖ To generate a compatible result file report

Set the ProSightPD FDR node's Generate tdReport parameter, to **True**. This node is available in the consensus workflow.

For more information, see [FDR node](#).

Use TDValidator Lite

To launch the TDValidator Lite application, click **TDValidator Lite** in the Proteome Discoverer Results view.

The TDValidator Lite application provides a comprehensive overview of the matching fragment ions for a given proteoform and spectrum. The TDValidator Lite window opens with the associated precursor and fragmentation spectra loaded as well as the proteoform from the results table.

Note If the raw files are not in the location specified in the study, TDValidator Lite will not work.

The TDValidator Lite application generates theoretical isotopic distributions based on the exact chemical formula of the proteoform and corresponding fragment ions. These isotopic distributions are then fit to the experimental data. Only ions relevant to the proteoform of interest are considered. As such, the matching fragment ions in TDValidator Lite can be slightly different than those found in ProSight searches from Xtract-detected fragment ions.

In addition to matching terminal ions of a proteoform, the TDValidator Lite application can be used to look at matching internal fragment ions. Because the number of possible internal fragment ions quickly scales for larger proteoforms, several additional features are present in the TDValidator Lite application to help understand result output. An FDR feature is now included that generates randomized sequences of the same length as the proteoform. A distribution of matches for the randomized sequences can then be produced to understand how many fragment ions, terminal or internal, match on average to a random protein sequence. Another key part of the TDValidator Lite application is the fitter score. In general, scores greater than 0.7 can be considered fairly robust, but scores all the way down to 0.5 can still be acceptable. When looking at internal fragment ions, a higher score can be necessary to limit false positives, particularly for very dense spectral data or large proteoforms (>30 kDa). Lastly, ppm differences between experimental and theoretical ions can be used to limit the scope of the error differences that are allowed. For example, if terminal ions are generally clustering around 0 ppm mass difference but internal fragment ions are much more widely spaced, limiting the allowable mass differences can reduce the false positive internal fragment ions.


Filter results

Refer to the *Proteome Discoverer User Guide* for information on filtering your results.

Visualize results with charts and graphs

This section lists the visualization formats available for the ProSightPD application.

Histograms




Abundances


Bar Charts | Pie Charts | Venn Diagrams

ce: Proteoforms : Abundances (Normalized) ✓


- Proteoforms : Abundances
- Proteoforms : Abundances (Normalized)
- Proteoforms : Abundances (Grouped)
- Proteoforms : Abundances (by Bio. Rep.)



Volcano Plot



Scatter Plot



Chromatogram Traces view

The Chromatogram Traces view displays the following:


- The chromatographic profiles of individual features or peaks that the application uses in calculating label-free and precursor ion quantification results
- The exact peak and feature traces that are used to derive quantification values
- The feature traces that are summed up from their corresponding peaks

In addition, you can separately display the single peaks that are subsumed in feature traces to provide an even more detailed view. Closely overlapping chromatograms (or traces) in the view indicate whether proteoforms eluted at the same time from the column and therefore are the same proteoform. They confirm the validity of the quantification results and the origin of the consensus features from the same proteoform. Conversely, peaks appearing far apart in the view indicate potential false positives.

The Chromatogram Traces view is available only when you perform label-free quantification or precursor ion quantification, but what it displays in the view depends on the .pdResult file page that is open.

Abundance is the area under the curve for each trace. In [Figure 29](#), the Consensus Features page of the table associated with the selected proteoform group shows the abundances displayed in the Chromatogram Traces view.

Figure 29. Chromatogram Traces view



Workflow templates and analysis templates

This chapter provides information about the templates for top-down analysis.

Table 13 lists the analysis and workflow templates for study creation.

Table 13. Analysis Templates, Processing Workflow Templates, and Consensus Workflow Templates (Sheet 1 of 3)

Template Name	Description
Analysis Templates	Includes both processing workflow and consensus workflow.
PSPD 1 percent FDR LFQ for HI HI data.PDANALYSIS	Includes the discovery proteomics search nodes and the nodes required to perform label-free proteoform quantitation on high resolution data.
PSPD 1 percent FDR LFQ for Med HI data.PDANALYSIS	Includes the nodes for discovery proteomics searches and label-free proteoform quantitation on medium high resolution data.
PSPD Discovery Proteomics with FDR and PFR Annotation.	Includes the same nodes as the PSPD Discovery Proteomics with FDR, with the addition of the PFR annotator node in the consensus workflow.
PSPD Discovery Proteomics with FDR for Med Hi MS3 data.	Supports analysis of “complex-down” data. The initial fragmentation event generates monomers which are targeted for MS3 fragmentation. The sequence information is derived from the MS3 data (ie PrSMs come from MS3) This templates supports unresolved precursor signal (ie monomer isotopes are unresolved).
PSPD Hi Hi TopDownStandard Analysis Template	Designed to be used with the provided Pierce TopDown Standard database and RAW file.
PSPD Med HI Analysis Template	Built for data with isotopically unresolved precursors. Covers the largest proteoform space including all annotated forms (PTMs, SNPs, and so on) in your database (annotated proteoform search), unknown truncations (subsequence search) and unknown modifications through the wide mass (1000 Da) tolerance annotated proteoform search. This search is the most time consuming. Designed for traditional high throughput LC/MS based workflows. Results are filtered to 1% FDR.
PSPD Discovery Proteomics with FDR.pdAnalysis	Identifies all proteoforms found in the user’s database, including known modifications. This search does not find unknown truncated proteoforms. This template is designed for quick analysis of traditional high throughput LC/MS based workflows. Results are filtered to 1% FDR.

Table 13. Analysis Templates, Processing Workflow Templates, and Consensus Workflow Templates (Sheet 2 of 3)

Template Name	Description
PSPD Comprehensive Discovery Proteomics with FDR.pdAnalysis	Covers the largest proteoform space including all annotated forms (PTMs, SNPs, and so on) in your database (annotated proteoform search), unknown truncations (subsequence search) and unknown modifications through the wide mass (1000 Da) tolerance annotated proteoform search. This search is the most time consuming. Designed for traditional high throughput LC/MS based workflows. Results are filtered to 1% FDR.
PSPD Unresolved Proteoform Discovery Proteomics with FDR.pdAnalysis	Intended for data with non-isotopically resolved precursors. Typically, this search template is used for lower resolution precursor scans or when analyzing larger or highly charged proteoforms. Results are filtered to 1% FDR.
PSPD Truncation Subsequence Search with FDR.pdAnalysis	Searches for truncated sequences, including those with any known PTMs or modifications.
PSPD Single Proteoform.pdAnalysis	Searches for a single user-defined proteoform, which does not require a database. This workflow is suitable for identifying one species in a complex LC/MS experiment or to simply confirm the identify of a single species from a mixture of infused proteoforms.
Processing Workflow Templates	
PSPD MS2 Precursor with MS3 Fragmentation HI HI	Supports analysis of “complex-down” data. The initial fragmentation event generates monomers, which are targeted for MS3 fragmentation. The sequence information is derived from the MS3 data (PrSMs come from MS3) This template supports isotopically resolved signal (monomer isotopes are resolved).
PSPD Discovery Proteomics with FDR for Med HI MS3 data data	Supports analysis of “complex-down” data. The initial fragmentation event generates monomers, which are targeted for MS3 fragmentation. The sequence information is derived from the MS3 data (PrSMs come from MS3) This templates supports unresolved precursor signal (monomer isotopes are unresolved).
PSPD PTCR Discovery Proteomics with FDR.pdProcessingWF	Identifies larger proteoforms from PTCR data.
PSPD Discovery Proteomics with FDR for HI HI data.pdProcessingWF	Performs searches for high-resolution data.
PSPD Discovery Proteomics with FDR for Unresolved Med HI data.pdProcessingWF	Performs proteoform identification on data containing unresolved precursors. Must be used with FDR nodes in consensus workflow.
PSPD LFQ for HI HI data.pdProcessingWF	Performs searches for high resolution data with LFQ analysis. Must be used with LFQ nodes in consensus workflow for LFQ results.

Table 13. Analysis Templates, Processing Workflow Templates, and Consensus Workflow Templates (Sheet 3 of 3)

Template Name	Description
PSPD LFQ for Med HI data.pdProcessingWF	Performs searches for medium high data with LFQ analysis. For label-free results, use with LFQ nodes in the consensus workflow.
PSPD Comprehensive Discovery Proteomics with FDR for HI HI data.pdProcessingWF	Performs both annotated and subsequence searches for high resolution data (resolved precursors), Must be used with FDR nodes in consensus workflow.
PSPD Single Proteoform Search.pdProcessingWF	Performs single user-defined proteoform searches.
PSPD Truncation Search with FDR.pdProcessingWF	Performs subsequence search only for high resolution data. Must be used with FDR nodes in consensus workflow.
Consensus Workflow Templates	
PSPD 1 percent FDR Consensus.pdConsensusWF	Returns results filtered to 1% FDR.
PSPD 1 percent FDR with PFR Annotation	Returns results filtered to 1%FDR includes PFR annotation
PSPD 1 percent FDR LFQ for HI HI data.pdconsensuswf	Returns LFQ data with 1% FDR.
PSPD No FDR Consensus.pdConsensusWF	Returns results without FDR control.

ProSightPD nodes

This chapter provides information about ProSightPD nodes that are specific to top-down analysis.

Contents

- [Feature Detector nodes](#)
- [Search nodes](#)
- [cRAWler nodes](#)
- [ProSightPD Consensus nodes](#)

Feature Detector nodes

The following topics describe the feature detection nodes:

- [Med Res. Feature Detector](#)
- [Hi Res. Feature Detector](#)

Med Res. Feature Detector

The Med Res. Feature Detector node uses the sliding window and the kDecon algorithms to perform spectral deconvolution and measure all of the deconvoluted features and their quantitation traces. The parameters used in this node are essential to the label-free quantitation workflows. Set the parameters according to the data being analyzed.

The Sliding Window Averaging Width RT parameter establishes the RT for the sliding window average. This value should match the full width at half max of a representative peak from the analyzed data. Thermo Fisher Scientific recommends measuring the peak width of an average feature in your data set. Do not choose the most abundant peak or one that is at noise level.

Table 14 lists the node parameters.

Table 14. (Sheet 1 of 3) ProSightPD Med Res. Feature Detector node

Parameter	Description
1. kDecon	
Precursor Maximum Mass	Specifies the maximum precursor mass to be considered by kDecon.
Precursor Minimum Mass	Specifies the minimum precursor mass to be considered by kDecon.
Number of Results	The maximum number of kDecon results to return for each averaged scan. Results are filtered by S/N prior to being returned.
Number of Iterations	kDecon can be run multiple times. After each iteration, charge state distributions from each mass are removed from the spectrum to limit false positives during subsequent iterations.
S/N Result Cutoff	The geometric mean signal-to-noise ratio S/N is calculated from the individual S/N of each charge state in the charge state distribution of a detected mass.
Peak Tolerance in PPM	Peak Tolerance PPM
2. Sliding Window Parameters	
Averaging Width RT	Specifies the retention time, or the width, of the sliding window, in minutes. Reducing this value improves time resolution but reduces execution speed and possibly sensitivity. Increasing this value increases execution speed but reduces time resolution and might increase sensitivity. You can achieve the best results when the window width is between one-quarter and twice the width of the characteristic peaks in the spectrum. For most usage, the optimum value might be half the width of the characteristic peaks. For example, if those peaks have a width of one minute, the optimum width would be 0.5 minutes.
Offset Type	Specifies the offset between successive sliding windows as a number of scans or as a percentage value. You can select from these two options: scan-based offset or percentage of averaging width RT offset.
Offset Scan(s)	This mode offsets each window from its predecessor by the user-specified number of scans. An offset of n means that each window begins n scans after the beginning of its predecessor.

Table 14. (Sheet 2 of 3) ProSightPD Med Res. Feature Detector node

Parameter	Description
Offset Percentage	This mode offsets each window from its predecessor by the user-specified percentage of the window width. An offset of 30% means that each window begins 30% after the beginning of its predecessor and overlaps the last 70%. An offset of 100% means that successive windows are adjacent with no overlap.
Merge Tolerance (ppm)	Determines how close two components in successive sliding windows must be in mass (ppm) for the application to identify them as a single component. A value that is too high might result in merging components that should remain separate. A value that is too low might result in false positives when components remain separate that should have merged.
Minimum Number of Charge States Observed	Minimum required number of observed charge states for data to be analyzed by KDecon. For top-down data, the recommended minimum is 3. Some very low-abundance species might only have 2 observable (greater than noise) peaks. However, setting this value to less than 3 will allow noise through and result in false positives.
Minimum Number of Sliding Window Detections	Specifies the minimum number of sliding window intervals that a component must appear in for the application to consider the component valid. A value that is too low might allow noise peaks to appear as false positives. A value that is too high might result in legitimate components being discarded. Set this parameter to a value large enough to exclude results that are implausibly narrow in retention time but small enough to include results of realistic duration. Values in the range of 3 through 8 generally give good results. A good approach is to use whichever is larger: 3 or the minimum number of windows that can fit into a peak.
Biggest Gap Type	Specifies the allowed separation between successive individual members of a merged component identified by the Sliding Window Deconvolution algorithm. You can select from these two options: Scans or Retention Time.
Biggest Gap Scans	The maximum allowed separation in scans between two successive individual members of a merged component identified by the Sliding Window Deconvolution algorithm. If the separation exceeds this value, the component will be divided into two components separated. As with the Merge Tolerance, too high a value can result in components that should have remained distinct being merged, while too low a value can result in false positives when components that should have been merged remain distinct. This parameter should be comparable too or slightly less than the expected separation in scans between distinct components with the same mass.

Table 14. (Sheet 3 of 3) ProSightPD Med Res. Feature Detector node

Parameter	Description
Biggest Gap Retention Time	The maximum allowed separation in retention time between two successive individual members of a merged component identified by the Sliding Window Deconvolution algorithm. If the separation exceeds this value, the component will be divided into two components separated by the observed gap in retention time. As with the Merge Tolerance, too high a value can result in components that should have remained distinct being merged, while too low a value can result in false positives when components that should have been merged remain distinct. This parameter should be comparable too or slightly less than the expected separation in retention time between distinct components with the same mass.
3. Feature Grouping	
Mass Tolerance	Mass tolerance used for calculating mass traces.
RT Threshold	Retention time limits used for calculating features.
4. Feature/Trace Connection	
Trace Mass Tolerance	Specifies the trace tolerance.
Number of Smoothing Points	Specifies the number of points to average. As the number increases, the data become smoother resulting in a loss of fine features. The recommended number of points is 3.
Time Range (min)	The allowable time difference between a trace and its linked feature. The recommended setting is 1 minute.
Trace Smoothing Type	Specifies the type of smoothing to perform: Gaussian, Moving Average, or None. The recommended setting is None. This should be used only in cases where unsmoothed data is not acceptable.
5. Multithreading Options	
CPU Usage	Amount of CPU to direct towards data processing.

Hi Res. Feature Detector

The Hi Res. Feature Detector node uses the sliding window and Xtract algorithms to perform spectral deconvolution and measure all of the deconvoluted features and their quantitation traces. The parameters used in this node are essential to the LFQ workflow. Set them according to the data you are analyzing.

One of the key parameters is the Sliding Window Averaging Width RT. This parameter establishes the RT for the sliding window average. Set this parameter to match the full width at half max of a representative peak from the analyzed data. It is recommended to measure the peak width of an average feature in your data set. Do not choose the most abundant peak or one that is at noise level.

The best parameter setting is based on your chromatography:

- If the window is too wide, then low intensity signals are averaged out of the data.
- If the window is too narrow, then you lose the advantages of signal averaging.

Table 15 lists the node parameters.

Table 15. ProSightPD Hi Res. Feature Detector parameters (Sheet 1 of 6)

Parameter	Description
1. Xtract	
S/N Threshold	Specifies the minimum signal-to-noise ratio for data that the node analyses.
Lowest m/z	Specifies the portion of the input spectrum that the Xtract algorithm processes. Min: Specifies the lowest end of the input spectrum.
Highest m/z	Specifies the portion of the input spectrum that the Xtract algorithm processes. Max: Specifies the highest end of the input spectrum.
Min Precursor Mass	The minimum precursor mass that the Xtract algorithm considers.
Max Precursor Mass	The maximum precursor mass that the Xtract algorithm considers.
Lowest Charge	Sets the low end of the allowable range for the number of charge states that must appear for a component to be recognized. The Xtract algorithm rejects potential components with fewer than the minimum or greater than the maximum number of charge states.
Highest Charge	Sets the high end of the allowable range.
Minimum Number of Detected Charge States	Minimum required number of observed charge states for data to be analyzed by the Xtract algorithm. For top-down data, the recommended minimum is 3. Some very low-abundance species might have only 2 observable (greater than noise) peaks. However, setting this value to less than 3 lets noise through and results in false positives.

Table 15. ProSightPD Hi Res. Feature Detector parameters (Sheet 2 of 6)

Parameter	Description
Relative Abundance Threshold	<p>Specifies a threshold below which the node filters out data for data reporting. This option sets a relative threshold as a percentage of the most abundant component in the spectrum. The most abundant peak in the deconvolved spectrum has a relative abundance of 100 percent, and all other peaks are calculated relative to that one. For example, if the highest peak has an absolute abundance of 1000, the relative abundance is 1 percent, and no peaks below an absolute abundance of 10 appear in the deconvolved spectrum.</p> <p>Range: 0–100</p> <p>Zero (0) displays all results; 100 displays only the most abundant component.</p>
Resolution at m/z 400	<p>Defines the resolution of the source spectrum at an m/z value of 400.</p>
Fit Factor	<p>Measures the quality of the match between a measured isotope pattern and an average distribution of the same mass. Enter a value between 0 and 100%</p> <ul style="list-style-type: none">• 0% requires a low fit only.• 100% means that the measured isotope profile is identical to the theoretical average isotope distribution.
Remainder Threshold (%)	<p>Specifies the height of the smaller overlapping isotopic cluster, as a percentage, with respect to the height of the most abundant isotopic cluster when the Xtract algorithm attempts to resolve overlapping isotopic clusters.</p> <p>For example, if one isotopic cluster in a spectrum has an abundance of 100, and you set the Remainder Threshold parameter to 30%, the Xtract algorithm ignores any overlapping clusters with an abundance less than 30.</p>
Min. Intensity	<p>Specifies a minimum intensity threshold to filter out possible background noise, including when you set the S/N Threshold parameter to zero.</p>
Expected Intensity Err.	<p>Specifies the permissible percentage of error allowed in calculating the ratio of the most abundant isotope to the next isotope higher in mass in the isotope series.</p>

Table 15. ProSightPD Hi Res. Feature Detector parameters (Sheet 3 of 6)

Parameter	Description
Consider Overlaps	<p>When selected (default), indicates that the Xtract algorithm is more tolerant of errors when the spectrum intensity is significantly higher than expected for the theoretical isotopic cluster.</p> <p>This option can lead to increased false positives; select it only when you expect overlapping isotopic clusters in a data set.</p>
2. Sliding Window Parameters	
Averaging Width RT	<p>Specifies the retention time, or the width, of the sliding window, in minutes.</p> <ul style="list-style-type: none"> Reducing this value improves time resolution but reduces execution speed and possibly sensitivity. Increasing this value increases execution speed but reduces time resolution and might increase sensitivity. <p>You can achieve the best results when the window width is between one-quarter and twice the width of the characteristic peaks in the spectrum. For most uses, the optimum value might be half the width of the characteristic peaks. For example, if those peaks have a width of one minute, the optimum width is 0.5 minutes.</p>
Offset Type	<p>Specifies the offset between successive sliding windows as a number of scans or as a percentage value. Select from:</p> <ul style="list-style-type: none"> Scan—Uses the Offset Scan setting Percent—Uses the Offset Percentage setting
Offset Scan	<p>Offsets each window from its predecessor by the user-specified number of scans. An offset of n means that each window begins n scans after the beginning of its predecessor.</p>
Offset Percentage	<p>Offsets each window from its predecessor by the user-specified percentage of the window width.</p> <ul style="list-style-type: none"> An offset of 30% causes each window to begin 30% after the beginning of its predecessor and overlap the last 70%. An offset of 100% causes successive windows to be adjacent with no overlap.

Table 15. ProSightPD Hi Res. Feature Detector parameters (Sheet 4 of 6)

Parameter	Description
Merge Tolerance	<p>Determines how close two components in successive sliding windows must be in mass for the application to identify them as a single component.</p> <ul style="list-style-type: none">• A value that is too high might result in merging components that should remain separate.• A value that is too low might result in false positives when components remain separate that should have merged.
Merge Mode	<p>Determines how spectra are merged:</p> <ul style="list-style-type: none">• Batch—Merges spectra as a batch, which minimizes the effect of scan to scan variation (Recommended).• Incremental—Merges spectra incrementally, one-by-one, making the result sensitive to scan-by-scan variations for long eluting species (Legacy).
Minimum Number of Charge States Observed	<p>Minimum required number of observed charge states for data to be analyzed by kDecon.</p> <p>For top-down data, the recommended minimum is 3. Some very low-abundance species might have only 2 observable (greater than noise) peaks. However, setting this value less than 3 lets noise through and result in false positives.</p>
Minimum Number of Sliding Window Detections	<p>Specifies the minimum number of sliding window intervals that a component must appear in for the application to consider the component valid.</p> <p>Set this parameter to a value large enough to exclude results that are implausibly narrow in retention time but small enough to include results of realistic duration. Values in the range of 3 through 8 generally give good results. A good approach is to use whichever is larger: 3 or the minimum number of windows that can fit into a peak.</p> <ul style="list-style-type: none">• A value that is too low might allow noise peaks to appear as false positives.• A value that is too high might result in legitimate components being discarded.

Table 15. ProSightPD Hi Res. Feature Detector parameters (Sheet 5 of 6)

Parameter	Description
Biggest Gap Type	Specifies the maximum allowed separation in retention time between two successive individual members of a merged component identified by the sliding windows algorithm, in minutes. <ul style="list-style-type: none"> Scans Retention time
Biggest Gap Scans	Specifies the greatest number of scans between subsequent feature detection for grouping similar features. <ul style="list-style-type: none"> If the separation exceeds this value, the application divides the candidate component into two merged components separated by a gap in retention time. A value that is too high for this parameter might result in merging components that should remain separate. A value that is too low might result in false positives when components remain separate that should have merged. <p>Set this parameter comparable to or slightly less than the expected separation in retention time between distinct components with the same mass.</p>
Biggest Gap Retention Time	Specifies the greatest amount of time between feature detection events for grouping similar features.
3. Feature /PrSM Connection	
Mass Tolerance	Delta mass (in Daltons) allowed for mapping PrSMs to features.
RT Threshold	Retention time threshold for mapping PrSMs to features: <ul style="list-style-type: none"> Minimum value = 0 Maximum value = (unchecked)
4. Feature/Trace Connection	
Trace Mass Tolerance	Specifies the trace tolerance. <ul style="list-style-type: none"> Minimum value = 0.01 Da 0.1 ppm Maximum value = (unchecked)
Number Of Smoothing Points	<ul style="list-style-type: none"> Minimum value = 0 Maximum value = (unchecked)
Time Range (min)	<ul style="list-style-type: none"> Minimum value = 0 Maximum value = (unchecked)

Table 15. ProSightPD Hi Res. Feature Detector parameters (Sheet 6 of 6)

Parameter	Description
Trace Smoothing Type	Determines the method to smooth the traces: <ul style="list-style-type: none">• Gauss• MovingMean (Default)• None
5. Multithreading Options	
CPU Usage	Determines the amount of CPU to direct toward processing: <ul style="list-style-type: none">• High—All available threads less 1 (Default)• Medium—50% of the available threads• Low—1 thread

Search nodes

This section describes the nodes that search the assigned database.

Note The ProSightPD 4.1 application includes the legacy search nodes from ProSightPD 3.0 so that you can reproduce previous results in the ProSightPD 4.1 application. These legacy nodes require a database in the .pscw file format generated in ProSightPC. You cannot directly convert .pscw to .psdb.

Use the source .xml file, which the .pscw file was generated from, to generate a .psdb-formatted database. To use a legacy search node, replace the 4.1 node with the equivalent legacy node and review all parameters to ensure agreement between your previous ProSightPD 3.0 workflows and your ProSightPD 4.1 workflows.

- [Annotated Proteoform Search node](#)
- [Single Proteoform Search node](#)
- [Subsequence Search node](#)
- [Tag to Annotated Proteoform Search node](#)
- [Tag to Subsequence Search](#)
- [ProSightPC PUF Writer node](#)

Annotated Proteoform Search node

The Annotated Proteoform Search node performs a ProSight-annotated proteoform search.

[Table 16](#) lists the parameters for the 4.1 version of the node. For information about the 3.0 version, refer to the user interface.

Table 16. ProSightPD 4.1 Annotated Proteoform Search node parameters (Sheet 1 of 3)

Parameter	Description
1. Input Data	
Input Database	Select the name of the database to search.
2. Annotated Proteoform (Absolute Mass) Search	
Precursor Mass Tolerance	Specifies the mass tolerance in Da or ppm used to find peak candidates. Range: 0.0–1 000 000 Da; 0.1–1000 ppm
Fragment Mass Tolerance	Specifies the mass tolerance in Da or ppm used for matching fragment peaks. Range: 0–2.0 Da; 1–50 ppm

Table 16. ProSightPD 4.1 Annotated Proteoform Search node parameters (Sheet 2 of 3)

Parameter	Description
Maximum PTMs per Isoform	Limits the number PTMs to add to any given proteoform. If more than this number of PTMs are present on an isoform, the proteoforms created are limited to the ones with combinations of PTMs having this many PTMS or fewer.
Maximum SNPs per Isoform	Limits the number SNPs to add to any given proteoform. If more than this number of SNPs are present on an isoform, the proteoforms created are limited to the ones with combinations of SNPs having this many SNPs or fewer.
Remove Disulfide Bonds	Thermo Fisher Scientific recommends setting this parameter to False. You can set this parameter to True if there are annotated disulfide bonds in the search database that you want to avoid searching, or when there is a known error in the disulfide bond annotation in the database. Errors in UniProt where disulfide bonds are placed on non-cysteine residues have been observed and can cause the search to crash.
Maximum Mass to Include PTMs	Maximum proteoform mass to search with PTMs.
N-Term Modifications to Include	Selected modifications are included on all intact proteoforms n-termini.
Delta M Mode	Determines whether the application conducts the search in Delta M (Δm) Mode. For more information on this mode, see Delta M Mode searches .
Calculate FDR	Specifies whether the FDR is calculated. You might see fewer results passing the FDR threshold when you use Delta M Mode.
Decoy Reps	The number of times a decoy search is run. A higher number of repetitions produces more stable decoy results but increases search times.
Maximum PrSMs per Precursor	Maximum number of PrSMs matched per precursor limits the number of results to the highest scoring hits.
Minimum Matched Fragments	Return only hits with this many or more matched fragments.
3. Static (Fixed) Modifications	
Static (Fixed) Modification	Specifies a static (also known as fixed) modification from a list of known modifications. Static modifications apply universally to every instance of the specified residues.

Table 16. ProSightPD 4.1 Annotated Proteoform Search node parameters (Sheet 3 of 3)

Parameter	Description
4. Multithreading Options	
CPU Usage	The amount of CPU to direct to data processing: <ul style="list-style-type: none"> • High—The total number of CPU cores less 1 • Medium—Half the number of available cores • Low—1 core
5. Fragmentation Method	
Override Fragmentation Method	Manually set the fragmentation method. Auto: allows the software to determine the correct fragmentation type (Default).
6. UVPD Method	
Default UVPD Method	Set the number of ion types to consider for the score of UVPD fragment ion spectra: <ul style="list-style-type: none"> • 4—recommended for maximizing IDs • 9—recommended for maximizing characterization

Single Proteoform Search node

The Single Proteoform Search node performs proteoform searches against a single proteoform entry. This mode also works for data without precursor scans (for example, direct infusion data). To run a search, the application requires a protein sequence and the RESID entry for any modified residues. The sequence format is described in [Table 17](#).

For an example of using this node in a workflow, see [Perform Single Proteoform searches in the processing workflow](#).

[Table 17](#) lists the parameters for the 4.1 version of the node. For information about the 3.0 version, refer to the user interface.

Table 17. ProSightPD 4.1 Single Proteoform Search parameters (Sheet 1 of 2)

Parameter	Description
1. Input Data	
RESID sequence	<p>Input a RESID sequence by inputting the single character abbreviation for each residue. You can add modifications to a residue by putting a RESID number in parentheses before the residue to modify.</p> <p>For example, a human Histone H3 proteoform with its N-terminal methionine removed modified by a K4 methylation and K36 trimethylation is input as:</p> <pre>ART(76)KQTARKSTGGKAPRKQLATKAARKSAPATGGV (74)KKPHRYRPGTVALREIRRYQKSTELLIRKLPFQRLVREI AQDFKTDLRFQSSAVMALQEACEAYLVGLFEDTNLCAIHA KRVTIMPKDIQLARRIRGERA</pre> <p>For the list of RESID modifications, go to: https://proteininformationresource.org/</p>
Include N-Term Acetylation	<p>Applies N-term acetylation to the RESID sequence:</p> <ul style="list-style-type: none"> • True • False
Proteoform Accession	Input an identifier that will appear in the pdResult file.
Proteoform Description	Input a description that will appear in the pdResult file.
2. Single Proteoform Search	
Precursor Mass Tolerance	<p>Specifies the mass tolerance in Da or ppm used to find peak candidates.</p> <p>Range: 0.0–1 000 000 Da; 0.1–1000 ppm</p>
Fragment Mass Tolerance	<p>Specifies the maximum mass difference for an observed fragment mass to be considered a match to a theoretical fragment mass:</p> <ul style="list-style-type: none"> • Minimum value = 1 ppm • Maximum value = 2 Da 50 ppm
Delta M Mode	Determines whether the application conducts the search in Delta M (Δm) Mode. For more information on this mode, see Delta M Mode searches .
Static (Fixed) Modification	Specifies a static (also known as fixed) modification from a list of known modifications. Static modifications apply universally to every instance of the specified residues.
Maximum PrSMs per Precursor	Maximum number of PrSMs matched per precursor limits the number of results to the highest scoring hits.

Table 17. ProSightPD 4.1 Single Proteoform Search parameters (Sheet 2 of 2)

Parameter	Description
Minimum Matched Fragments	Return only hits with this many or more matched fragments.
3. Static (Fixed) Modifications	
Static (Fixed) Modification	Specifies a static (also known as fixed) modification from a list of known modifications. Static modifications apply universally to every instance of the specified residues.
4. Multithreading Options	
CPU Usage	The amount of CPU to direct to data processing: <ul style="list-style-type: none"> • High—The total number of CPU cores less 1 • Medium—Half the number of available cores • Low—1 core
5. Fragmentation Method	
Override Fragmentation Method	Manually set the fragmentation method. Auto: allows the software to determine the correct fragmentation type (Default).
6. UVPD Method	
Default UVPD Method	Set the number of ion types to consider for the score of UVPD fragment ion spectra: <ul style="list-style-type: none"> • 4—recommended for maximizing IDs • 9—recommended for maximizing characterization
7. Output Data	
Export Results to ProSightPC	Determines is a PUF file is created with the search results.

Subsequence Search node

The Subsequence Search node performs a subsequence search. Increasing the precursor tolerance results in longer run times.

To search for modified proteoforms, select the Include Modified Forms check box in the Predefined Search dialog box for biomarkers.

[Table 18](#) lists the parameters for the 4.1 version of the node. For information about the 3.0 version, refer to the user interface.

Table 18. ProSightPD 4.1 Subsequence Search parameters (Sheet 1 of 2)

Parameter	Description
1. Input Data	
Input Database	Specifies the name of the database to search.
2. Subsequence (Biomarker) Search	
Subsequence Precursor Mass Tolerance	Specifies the maximum mass difference between the subsequence mass and the observed proteoform mass for a proteoform to be included in the search. Setting this to a large tolerance can potentially increase search time with large databases. Minimum value = 0.01 Da (10.1 ppm) Maximum value = 1 000 000 Da (1000 ppm)
Fragment Mass Tolerance	Specifies the mass tolerance in Da or ppm used for matching fragment peaks. Range: 0–2.0 Da; 1–50 ppm
Maximum PTMs per Isoform	Limits the number PTMs to add to any given proteoform. If more than this number of PTMs are present on an isoform, the proteoforms created are limited to the ones with combinations of PTMs having this many PTMS or fewer.
Maximum SNPs per Isoform	Limits the number SNPs to add to any given proteoform. If more than this number of SNPs are present on an isoform, the proteoforms created are limited to the ones with combinations of SNPs having this many SNPs or fewer.
Remove Disulfide Bonds	Thermo Fisher Scientific recommends setting this parameter to False. You can set this parameter to True if there are annotated disulfide bonds in the search database that you want to avoid searching, or when there is a known error in the disulfide bond annotation in the database. Errors in UniProt where disulfide bonds are placed on non-cysteine residues have been observed and can cause the search to crash.
Maximum Mass to Include PTMs	Maximum proteoform mass to search with PTMs.
N-Term Modifications to Include	Selected modifications are included on all intact proteoforms n-termini.
Calculate FDR	Specifies whether the FDR is calculated. You might see fewer results passing the FDR threshold if you use Delta M Mode.

Table 18. ProSightPD 4.1 Subsequence Search parameters (Sheet 2 of 2)

Parameter	Description
Decoy Reps	The number of times a decoy search is run. A higher number of repetitions produces more stable decoy results but increases search times.
Maximum PrSMs per Precursor	Maximum number of PrSMs matched per precursor limits the number of results to the highest scoring hits.
Minimum Matched Fragments	Return only hits with this many or more matched fragments.
3. Static (Fixed) Modifications	
Static (Fixed) Modification	Specifies a static (also known as fixed) modification from a list of known modifications. Static modifications apply universally to every instance of the specified residues.
4. Multithreading Options	
CPU Usage	The amount of CPU to direct to data processing: <ul style="list-style-type: none"> • High—The total number of CPU cores less 1 • Medium—Half the number of available cores • Low—1 core
5. Fragmentation Method	
Override Fragmentation Method	Manually set the fragmentation method. Auto: allows the software to determine the correct fragmentation type. (Default)
6. UVPD Method	
Default UVPD Method	Set the number of ion types to consider for the score of UVPD fragment ion spectra: <ul style="list-style-type: none"> • 4—recommended for maximizing IDs • 9—recommended for maximizing characterization

Tag to Annotated Proteoform Search node

The Tag to Annotated Proteoform Search node performs a sequence tag search first and then sends all the results from the sequence tag search to an annotated proteoforms search. The sequence tag search reduces the amount of spectra searched by the annotated proteoform search step and reduces the overall search time for large data.

[Table 19](#) lists the parameters for the 4.1 version of the node. For information about the 3.0 version, refer to the user interface.

Table 19. ProSightPD 4.1 Tag to Annotate Proteoform Search parameters (Sheet 1 of 2)

Parameter	Description
1. Input Data	
Input Database	Specifies the name of the database to search.
2a. Sequence Tag	
Minimum Tag Score	Minimum probability that the probable tag is not found by chance.
Amino Acid Gap	Largest number allowed of sequential amino acids not having a fragment between them.
Max Tags	Greatest number of probable tags considered for each search per target.
Minimum Fragments Used	Minimum number of fragments required for each sequence tag.
Remove Water Losses	When set to True, water loss is not considered when searching sequence tags.
Consider Isobaric Dipeptide Gaps	When set to True, for example, both GA and Q are considered for gaps of 128.0 Da.
Compiler Tolerance	The permissible error, measured in ppm between two fragment ion masses that are still considered as matching an amino acid
2b. Annotated Proteoform (Absolute Mass) Search	
Fragment Mass Tolerance	Specifies the mass tolerance in Da or ppm used for matching fragment peaks. Range: 0–2.0 Da; 1–50 ppm
Maximum PTMs per Isoform	Limits the number PTMs to add to any given proteoform. If more than this number of PTMs are present on an isoform, the proteoforms created are limited to the ones with combinations of PTMs having this many PTMS or fewer.
Maximum SNPs per Isoform	Limits the number SNPs to add to any given proteoform. If more than this number of SNPs are present on an isoform, the proteoforms created are limited to the ones with combinations of SNPs having this many SNPs or fewer.
Remove Disulfide Bonds	Thermo Fisher Scientific recommends setting this parameter to False. You can set this parameter to True if there are annotated disulfide bonds in the search database that you want to avoid searching, or when there is a known error in the disulfide bond annotation in the database. Errors in UniProt where disulfide bonds are placed on non-cysteine residues have been observed and can cause the search to crash.

Table 19. ProSightPD 4.1 Tag to Annotate Proteoform Search parameters (Sheet 2 of 2)

Parameter	Description
N-Term Modifications to Include	Selected modifications are included on all intact proteoforms n-termini.
Delta M Mode	Determines whether the application conducts the search in Delta M (Δm) Mode. For more information on this mode, see Delta M Mode searches .
Calculate FDR	Specifies whether the FDR is calculated. You might see fewer results passing the FDR threshold if you use Delta M Mode.
Decoy Reps	The number of times a decoy search is run. A higher number of repetitions produces more stable decoy results but increases search times.
Maximum PrSMs per Precursor	Maximum number of PrSMs matched per precursor limits the number of results to the highest scoring hits.
Minimum Matched Fragments	Return only hits with this many or more matched fragments.
3. Static (Fixed) Modifications	
Static (Fixed) Modification	Specifies a static (also known as fixed) modification from a list of known modifications. Static modifications apply universally to every instance of the specified residues.
4. Multithreading Options	
CPU Usage	The amount of CPU to direct to data processing: <ul style="list-style-type: none"> • High—The total number of CPU cores less 1 • Medium—Half the number of available cores • Low—1 core
5. Fragmentation Method	
Override Fragmentation Method	Manually set the fragmentation method. Auto: allows the software to determine the correct fragmentation type. (Default)
6. UVPD Method	
Default UVPD Method	Set the number of ion types to consider for the score of UVPD fragment ion spectra: <ul style="list-style-type: none"> • 4—recommended for maximizing IDs • 9—recommended for maximizing characterization

Tag to Subsequence Search

The Tag to Subsequence Search node first performs a tag search and then performs a sequence tag search on the output of the tag search.

Table 20 lists the parameters for the 4.1 version of the node. For information about the 3.0 version, refer to the user interface.

Table 20. ProSightPD 4.1 Tag to Subsequence Search node parameters (Sheet 1 of 3)

Parameter	Description
1. Input Data	
Input Database	Select the search database.
2a. Subsequence Tag	
Minimum Tag Score	Minimum probability that the probable tag is not found by chance.
Amino Acid Gap	Largest number allowed of sequential amino acids not having a fragment between them.
Max Tags	Greatest number of probable tags considered for each search per target.
Minimum Fragments Used	Minimum number of fragments required for each sequence tag.
Remove Water Losses	If set to True, water loss is not considered when searching sequence tags.
Consider Isobaric Dipeptide Gaps	If set to True, for example, both GA and Q are considered for gaps of 128.0 Da.
Compiler Tolerance	The permissible error, measured in ppm between two fragment ion masses that are still considered as matching an amino acid.
2b. Subsequence (Biomarker) Search	
Subsequence Precursor Mass Tolerance	Specifies the tolerance when comparing subsequence mass to observed precursors. A large tolerance setting might increase search time with large databases. <ul style="list-style-type: none"> • Minimum value = 0.01 Da 10.1 ppm • Maximum value = 1 000 000 Da 1000 ppm
Fragment Mass Tolerance	Specifies the mass tolerance in Da or ppm used for matching fragment peaks. Range: 0–2.0 Da; 1–50 ppm

Table 20. ProSightPD 4.1 Tag to Subsequence Search node parameters (Sheet 2 of 3)

Parameter	Description
Maximum PTMs per Isoform	Limits the number PTMs to add to any given proteoform. If more than this number of PTMs are present on an isoform, the proteoforms created are limited to the ones with combinations of PTMs having this many PTMS or fewer.
Maximum SNPs per Isoform	Limits the number SNPs to add to any given proteoform. If more than this number of SNPs are present on an isoform, the proteoforms created are limited to the ones with combinations of SNPs having this many SNPs or fewer.
Maximum Mass to Include PTMs	Maximum proteoform mass to search with PTMs.
Remove Disulfide Bonds	Thermo Fisher Scientific recommends setting this parameter to False. You can set this parameter to True if there are annotated disulfide bonds in the search database that you want to avoid searching, or when there is a known error in the disulfide bond annotation in the database. Errors in UniProt where disulfide bonds are placed on non-cysteine residues have been observed and can cause the search to crash.
N-Term Modifications to Include	Selected modifications are included on all intact proteoforms n-termini.
Delta M Mode	Determines whether the application conducts the search in Delta M (Δm) Mode. For more information on this mode, see Delta M Mode searches .
Calculate FDR	Specifies whether the FDR is calculated. You might see fewer results passing the FDR threshold if you use Delta M Mode.
Decoy Reps	The number of times a decoy search is run. A higher number of repetitions produces more stable decoy results but increases search times.
Maximum PrSMs per Precursor	Maximum number of PrSMs matched per precursor limits the number of results to the highest scoring hits.
Minimum Matched Fragments	Return only hits with this many or more matched fragments.
3. Static (Fixed) Modifications	
Static (Fixed) Modification	Specifies a static (also known as fixed) modification from a list of known modifications. Static modifications apply universally to every instance of the specified residues.

Table 20. ProSightPD 4.1 Tag to Subsequence Search node parameters (Sheet 3 of 3)

Parameter	Description
4. Multithreading Options	
CPU Usage	The amount of CPU to direct to data processing: <ul style="list-style-type: none"> • High—The total number of CPU cores less 1 • Medium—Half the number of available cores • Low—1 core
5. Fragmentation Method	
Override Fragmentation Method	Manually set the fragmentation method. Auto: allows the software to determine the correct fragmentation type. (Default)
6. UVPD Method	
Default UVPD Method	Set the number of ion types to consider for the score of UVPD fragment ion spectra: <ul style="list-style-type: none"> • 4—recommended for maximizing IDs • 9—recommended for maximizing characterization

ProSightPC PUF Writer node

The ProSightPC PUF Writer node is a legacy search node. For more information, refer to the node's parameter listing in the user interface.

cRAWler nodes

This section describes the cRAWler nodes:

- [APD cRAWler node](#)
- [High/High cRAWler node](#)
- [Med/High cRAWler node](#)

APD cRAWler node

The APD cRAWler node processes PTCR precursor data with the APD algorithm and uses the Xtract algorithm to deconvolute fragmentation spectra.

[Table 21](#) lists the node parameters.

Table 21. ProSightPD APD cRAWler node parameters (Sheet 1 of 2)

Parameter	Description
1. Chromatographic Parameters	
Max RT Difference [min]	The maximum time elapsed between MS/MS scans to be grouped together into a target.
Precursor m/z Tolerance	The maximum difference in precursor isolation m/z allowed for grouping of subsequent MS/MS scans.
Grouping	Defines the grouping method: <ul style="list-style-type: none"> • Standard—Groups MS/MS scans within retention time and m/z tolerances and groups the corresponding MS scans • PTCR—Groups a single fragmentation scan with its corresponding PTCR precursor scan
2. Precursor APD Parameters	
Use Manual Precursor Mass	Selects from the following: <ul style="list-style-type: none"> • Never: Always uses APD detected masses • Always: Only uses the manually input Precursor Mass parameter for every target mass • If not detected: When no mass is detected by APD the manually input Precursor Mass parameter is used for the target
Precursor Mass (Da) Override	User-defined precursor mass in Daltons.
Use Instrument Determined Precursor Mass	Allows the instrument to determine the precursor mass.
3. Fragmentation Xtract Parameters	
S/N Threshold	The minimum signal to noise threshold for a peak that the Xtract algorithm processes.
Lowest m/z	The lowest m/z considered by the Xtract algorithm.
Highest m/z	The highest m/z considered by the Xtract algorithm.
Lowest Charge	The lowest charge considered by the Xtract algorithm.
Highest Charge	The highest charge considered by the Xtract algorithm.
Fit Factor	The minimum pattern accuracy for an isotopic distribution to be considered a valid Xtract mass.
Overlapping Remainder	Minimum relative remaining intensity of an isotope pattern overlapped by another pattern to be considered.

Table 21. ProSightPD APD cRAWler node parameters (Sheet 2 of 2)

Parameter	Description
Resolution at m/z 400	Instrumental Resolution at m/z 400.
4. Multithreading Options	
CPU Usage	Determines the amount of CPU to direct toward processing: <ul style="list-style-type: none"> • High—All available threads less 1 (Default) • Medium—50% of the available threads • Low—1 thread
5. Data Reduction Parameters	
Enable Noise Reduction	Makes the following noise reduction parameters available.
Relative Multiple Precursor Threshold	Only consider additional precursors greater than the set value.
Absolute Minimum Fragment Intensity	Only considered fragments with intensity greater than the set value.
Window Size Below	Bin size below m/z 2000 for considering fragments.
Window Size Above	Bin size above m/z 2000 for considering fragments.
Fragments per window size above	Number of fragments per bin above m/z 2000 to consider.
Fragments per window size below	Number of fragments per bin below m/z 2000 to consider.
Minimum Fragment Count	Minimum number of fragments for a target to be searched.
Number of Precursors Per Experiment	Maximum number of precursors for a given experiment.

High/High cRAWler node

The High/High cRAWler node uses the Xtract algorithm to deconvolute the spectra from High/High data. The node then groups the deconvoluted spectra based on your parameter settings.

Table 22 lists the node parameters.

Table 22. ProSightPD High/High cRAWler node parameters (Sheet 1 of 3)

Parameter	Description
1. Chromatographic Parameters	
Max RT Difference [min]	The maximum time elapsed between MS/MS scans to be grouped together into a target.
Precursor m/z Tolerance	The maximum difference in precursor isolation m/z allowed for grouping of subsequent MS/MS scans.
Grouping	Determines the grouping method: <ul style="list-style-type: none"> • Standard—Groups MS/MS scans within retention time and m/z tolerances and groups the corresponding MS scans • PTCR—Groups a single fragmentation scan with its corresponding PTCR precursor scan
2. Precursor Xtract Parameters	
S/N Threshold	The minimum signal to noise threshold for a peak that the Xtract algorithm processes.
Lowest m/z	The lowest m/z considered by the Xtract algorithm.
Highest m/z	The highest m/z considered by the Xtract algorithm.
Lowest Charge	The lowest charge considered by the Xtract algorithm.
Highest Charge	Highest Charge considered by the Xtract algorithm.
Fit Factor	The minimum pattern accuracy for an isotopic distribution to be considered a valid Xtract mass.
Overlapping Remainder	Minimum relative remaining intensity of an isotope pattern overlapped by another pattern to be considered.
Resolution at m/z 400	Instrumental Resolution at m/z 400
Use Manual Precursor Mass	Selects from the following: <ul style="list-style-type: none"> • Never—Always uses the Xtract algorithm detected masses (Default) • Always—Only uses the manually input Precursor Mass parameter for every target mass • If not detected—When no mass is detected by the Xtract algorithm, the manually input Precursor Mass parameter is used for the target
Precursor Mass (Da) Override	User-defined precursor mass in Daltons.

Table 22. ProSightPD High/High cRAWler node parameters (Sheet 2 of 3)

Parameter	Description
Use Instrument Determined Precursor Mass	Sets precursor mass determined by the instrument.
3. Fragmentation Xtract Parameters	
S/N Threshold	The minimum signal-to-noise threshold for a peak that the Xtract algorithm processes.
Lowest m/z	The lowest m/z considered by the Xtract algorithm.
Highest m/z	The highest m/z considered by the Xtract algorithm.
Lowest Charge	The lowest charge considered by the Xtract algorithm.
Highest Charge	The highest charge considered by the Xtract algorithm.
Fit Factor	The minimum pattern accuracy for an isotopic distribution to be considered a valid Xtract mass.
Overlapping Remainder	Minimum relative remaining intensity of an isotope pattern overlapped by another pattern to be considered.
Resolution at m/z 400	Instrumental Resolution at m/z 400.
4. Multithreading Options	
CPU Usage	Determines the amount of CPU to direct toward processing: <ul style="list-style-type: none"> • High—All available threads less 1 (Default) • Medium—50% of the available threads • Low—1 thread
5. Data Reduction Parameters	
Enable Noise Reduction	Makes the following noise reduction parameters available.
Relative Multiple Precursor Threshold	Only consider additional precursors greater than the set value.
Absolute Minimum Fragment Intensity	Only considered fragments with intensity greater than the set value.
Window Size Below	Bin size below m/z 2000 for considering fragments
Window Size Above	Bin size above m/z 2000 for considering fragments.
Fragments per window size above	Number of fragments per bin above m/z 2000 to consider.
Fragments per window size below	Number of fragments per bin below m/z 2000 to consider.

Table 22. ProSightPD High/High cRAWler node parameters (Sheet 3 of 3)

Parameter	Description
Minimum Fragment Count	Minimum number of fragments for a target to be searched.
Number of Precursors Per Experiment	Maximum number of precursors for a given experiment.

Med/High cRAWler node

The Med/High cRAWler node provides mass assignments for precursor species in medium-high data with unresolved isotope distributions. The Med/High cRAWler node uses an iterative charge-state assignment approach to examine the most abundant m/z species in a spectrum and find real charge-state distributions. This approach is suitable for mass determination of larger mass precursor species.

Table 23 lists the node parameters.

Table 23. ProSightPD Med/High cRAWler Node parameters (Sheet 1 of 3)

Parameter	Description
1. Chromatographic Parameters	
Max. RT Difference[min]	The maximum time elapsed between MS/MS scans to be grouped together into a target.
Precursor m/z Tolerance	Specifies the maximum difference in precursor isolation m/z allowed for grouping of subsequent MS/MS scans.
Grouper	The grouping method: <ul style="list-style-type: none"> • Standard—Groups MS/MS scans within retention time and m/z tolerances and groups the corresponding MS scans. • PTCR—Groups a single fragmentation scan with its corresponding PTCR precursor scan.
2. Precursor kDecon Parameters	
Max Precursor Mass	Specifies the maximum precursor mass to consider.
Min Precursor Mass	Specifies the minimum precursor mass to consider.
Number Of Results	The maximum number of mass results to return for each averaged scan. The node filters the results by S/N before returning them.
Number of Iterations	The number of times to run the node. After each iteration, the node removes charge-state distributions from each mass contained in the spectrum to limit false positives during subsequent iterations.

Table 23. ProSightPD Med/High cRAWler Node parameters (Sheet 2 of 3)

Parameter	Description
S/N Result Cutoff	The geometric mean signal to noise ratio (S/N) is calculated from the individual S/N of each charge state in the charge state distribution of a detected mass.
Full Noise Reduction	Specifies whether to use a QuickHull algorithm to reduce noise. When set to False, the node removes peaks that are not above an exponential moving average.
Use Manual Precursor Mass	Selects from the following: <ul style="list-style-type: none"> • Never—Always uses algorithm detected masses (Default) • Always—Uses only the manually input Precursor Mass parameter for every target mass • If not detected—When the Xtract algorithm does not detect a mass, the manually input Precursor Mass parameter is used for the target
Precursor Mass (Da) Override	User-defined precursor mass in Daltons.
Use Instrument Determined Precursor Mass	Sets precursor mass determined by the instrument.
Charge State Intensity Type	The intensity returned for each mass is the sum of the detected charge states. If this parameter is set to Charge Normalized, the intensity of each charge state is divided by its charge.
Sum Centroid Intensities	If set to True, all intensity values that are part of the centroided peak are summed together. If set to False, only the most abundant m/z value is used.
3. Fragmentation Xtract Parameters	
S/N Threshold	The minimum S/N value for processing a peak.
Lowest m/z	Lowest m/z value to consider.
Highest m/z	Highest m/z value to consider.
Lowest Charge	The lowest charge considered by the Xtract algorithm.
Highest Charge	The highest charge considered by the Xtract algorithm.
Fit Factor	The minimum pattern accuracy for an isotopic distribution to be considered a valid Xtract mass.
Overlapping Remainder	Minimum relative remaining intensity of an isotope pattern overlapped by another pattern to be considered.
Resolution at m/z 400	Instrumental Resolution at m/z 400

Table 23. ProSightPD Med/High cRAWler Node parameters (Sheet 3 of 3)

Parameter	Description
4. Multithreading Options	
CPU Usage	Determines the amount of CPU to direct toward processing: <ul style="list-style-type: none"> • High—All available threads less 1 • Medium—50% of the available threads • Low—1 thread
5. Data Reduction Parameters	
Enable Noise Reduction	Makes the following noise reduction parameters available.
Relative Multiple Precursor Threshold	Consider only additional precursors greater than the set value.
Absolute Minimum Fragment Intensity	Consider only fragments with intensity greater than the set value.
Window Size Below	Bin size below m/z 2000 for considering fragments.
Window Size Above	Bin size above m/z 2000 for considering fragments.
Fragments per window size above	Number of fragments per bin above m/z 2000 to consider.
Fragments per window size below	Number of fragments per bin below m/z 2000 to consider.
Minimum Fragment Count	Minimum number of fragments for a target to be searched.
Number of Precursors Per Experiment	Maximum number of precursors for a given experiment.

ProSightPD Consensus nodes

This section describes all the top-down consensus nodes:

- [PrSM Grouper node](#)
- [Protein Grouper node](#)
- [Proteoform Validator node](#)
- [Proteoform and Protein Filter node](#)
- [Feature Mapper node](#)
- [Quantifier node](#)

- [FDR node](#)
- [PFR Annotator node](#)

PrSM Grouper node

The Proteoform Spectral Match Grouper node groups redundant PrSMs into proteoforms. In the consensus workflow, this node follows the MSF files node.

[Table 24](#) lists the advanced node parameters. The node does not have regular parameters.

Table 24. ProSightPD PrSM Grouper advanced node parameters

Parameter	Description
1. Proteoform Filters	
C-Score Threshold:	<p>Medium PrSM C-Score Threshold: PrSMs with C-Scores greater than or equal to this value, but less than the High PrSM C-Score Threshold are grouped together into medium confidence proteoforms.</p> <p>A proteoform with at least one PRSM of medium confidence is considered medium confidence.</p> <p>Default: 3</p>
High PrSM C-Score Threshold	<p>PrSMs with C-Scores greater than or equal to this value are grouped together into high confidence proteoforms.</p> <p>A proteoform with at least one PrSM of high confidence is considered high confidence.</p> <p>Default: 40</p>
2. Grouping by Tolerance	
Should group by tolerance	<p>If set to True, software groups proteoforms by mass in addition to sequence and modification.</p> <p>If set to False, software groups proteoforms only on their sequence and modifications. (Default)</p>
Group Detected Mass Tolerance	Specifies the grouping mass tolerance.

Protein Grouper node

The Protein Grouper node groups proteoforms and protein isoforms into their parent proteins. The node includes results thresholding and grouping options.

No parameters.

FDR node

You can use this node's parameter setting cutoff to include search results in the pdResult report.

The FDR calculations are implemented as described in the 2019 Mol Cell Proteomics article, “Accurate Estimation of Context–Dependent False Discovery Rates in Top–Down Proteomics,” by LeDuc et. al. For more information, go to <https://pubmed.ncbi.nlm.nih.gov/30647073>.

Table 25 lists the node parameters.

Table 25. ProSightPD FDR node parameters

Parameter	Description
1. FDR	
FDR cutoff	Enter a number that search results must exceed to be included in the pdResult report.
2. tdReport	
Generate tdReport	Generate a tdReport from the search results. The TD Viewer application can open the tdReport.

Proteoform Validator node

The Proteoform Validator node calculates proteoform confidences based on associated PrSMs.

No parameters.

Proteoform and Protein Filter node

The Proteoform and Protein Filter node filters proteoforms based on confidence.

Table 26 lists the node parameters.

Table 26. ProSightPD Proteoform and Protein Filter Node

Parameter	Description
1. Proteoform Filters	
Proteoform Confidence at Least	<p>Proteoforms with a lower confidence than specified are excluded from the final result:</p> <ul style="list-style-type: none"> • High (Default) • Medium • Low <p>You can define confidence levels in the configuration panel of the Administration page.</p>
Keep Lower Confidence PrSMs	<p>If set to True, keep PrSMs with lower confidence than specified at the proteoform level.</p> <p>If set to False, do not keep them. (Default)</p>

Feature Mapper node

The Feature Mapper node maps quantification consensus features onto proteoforms within the tolerances specified by the node.

Table 27 lists the node parameters.

Table 27. ProSightPD Feature Mapper node parameters

Parameter	Description
1. Consensus Feature	
Feature Presence in File Threshold	<p>Species minimum percentage of files a species must be quantified to return a quantitation result.</p> <p>For example, if the threshold is set to 30%, then a proteoform must be detected in at least 30% of the files analyzed to report a quantitation value.</p>
2. Consensus Features Group/Proteoform Connection	
Feature Group Proteoform Connection Time Tolerance (Minutes)	The maximum allowable separation in time for two feature groups to be grouped together.
Feature Group Proteoform Connection Mass Tolerance (Da or ppm)	The maximum allowable mass difference for two feature groups to be grouped together.

Quantifier node

The Quantifier node applies normalization and statistics to the resulting quantitative information. To use the volcano plots and heat maps, this node must be present in the consensus workflow.

Table 28 lists the node parameters.

Table 28. ProSightPD Quantifier node parameters

Parameter	Description
1. Settings	
Normalize	If set to True, normalization is applied to the proteoform level quantitation values.
Use Shared Results	If set to True, proteoforms can share quantification values. Due to minor differences in modification states (for example, positional isomers) frequently similar proteoforms coelute and therefore share elution profiles. If non-similar, non-coeluting species utilize shared quantitation values, Thermo Fisher Scientific recommends lowering the feature mapping tolerances or setting Used Shared Results to False.

PFR Annotator node

The PFR Annotator node assigns PFR Accessions to proteoforms that pass the FDR threshold in the ProSightPD FDR node. Only proteoforms associated with a protein entry and a Taxon Id will be annotated. Additionally, if there are multiple Taxon Ids associated with the search result, PFR Accessions will not be assigned. The PFR Annotator node adds a column proteoform results table, which includes a proteoform reference number (PFR). This lets you track proteoforms across different experiments.

No parameters.

Interpreting results

This chapter provides information about scoring systems shown in the Sequence Gazer for ProSightPD results.

Contents

- [P-Score](#)
- [C-Score](#)
- [Expectation Value \(E-Value\)](#)

P-Score

A P-Score is the probability of obtaining at least as good a match between the observed fragment list and a sequence as by chance. It is a measure of confidence in the validity of a match. A low P-Score means that the probability of obtaining at least this many fragments matching a sequence is low, so it is unlikely that random chance is the cause of the association.

The application calculates a P-Score as follows:

$$p(n) = \sum_{i=0}^{n-1} \frac{e^{-xf}(xf)^i}{i!}$$

where:

- n is the number of matching fragments.
- x is the probability of an observed fragment ion matching a random theoretical fragment ion by chance.
- f is the total number of fragment ions observed.

Matching more fragment ions leads to a lower, more confident P-Score. Alternatively, increasing the numbers of total detected fragment ions without a corresponding increase in matched fragment ions negatively affects the confidence of the result. To see the original score derivation, refer to:

<https://www.nature.com/articles/nbt1001-952>

C-Score

The C-Score measures the level of characterization of a proteoform in relation to the others in the database. This score, described in LeDuc et al.¹, uses a Bayesian approach that assigns a likelihood to each candidate proteoform based on the observed MS data.

A C-Score of 3 indicates that there are two proteoforms in the database that equally explain the observed data. A score of 40 or higher is considered strong evidence of a unique characterization.

Expectation Value (E-Value)

The expectation value (E-Value) is the number of sequences in a database that are expected to have P-Scores equal to or better than what was observed simply by chance. Low E-Values represent better matches (less likely to be false positives) than high E-Values. Because the P-Score represents the probability of the n out of f fragments matching by chance, and if it is assumed that all sequences in the database are independent, the E-Value of a sequence-fragment set association is simply the association's p value times the number of sequences in the database.

If N is the number of proteoforms considered during a search, the E-Value currently reported by the software is:

$$e = N \times p(n)$$

¹ LeDuc, R.D.; Fellers, R.T.; Early, B.P.; Greer, J.B.; Thomas, P.M.; Kelleher, N.L. The C-score: a Bayesian framework to sharply improve proteoform scoring in high-throughput top-down proteomics. *J Proteome Res.* 2014 Jul 3, 13 (7):3231-40.